

# Curera<sup>®</sup> Ear Pillow

For support and relief of sensitive ears



Ear pillow with cotton cover/hygiene cover

# Curera® Ear Pillow

A pillow designed to support and relieve sensitive ears, an area that is particularly prone to pressure ulcers when a user lies still in bed.

The pillow supports the face and the oval-shaped hole totally relieves pressure off the ear.

The ear pillow comes with a cotton cover and a hygiene cover. Accordingly, it ought to be used as a product intended for a single user.



The hygiene cover does not allow liquid or air to penetrate, but is vapour permeable (vapour is wicked away). Make sure that the patient is positioned correctly so that airways are unobstructed.



Read the user manual prior to using the product.

Care of Sweden and Curera are registered trademarks of Care of Sweden AB. © Care of Sweden AB, 2018. All rights reserved. Reproduction prohibited. Last updated: 09/11/2018.

CLEANING	FACTS
<p><b>Hygiene cover</b> The cover can be wiped clean with detergent and/or disinfectant. Machine wash at max. 95°C.</p> <p><b>Cotton cover</b> Machine wash at max. 70°C.</p>	<ul style="list-style-type: none"><li>• <b>Filling: 50% HR foam/50% polyester fibre</b></li><li>• <b>Hygiene cover: Polyurethane/polyester</b></li><li>• <b>Warranty: 1 year</b></li><li>• <b>Size: 32x32 cm</b></li></ul>