

Mattress Covers

A stack of several mattress covers in shades of blue and grey. Each cover has a white strap with a buckle, likely for securing the cover to the mattress. The covers are neatly folded and stacked on top of each other.

Meets health care service requirements
for hygiene, quality and safety

Hygiene covers (full cover or top part only)



Olivia^(1,2)

- Colour: grey.
- Zip on three sides.
- Composition: 55 % Polyester/45 % Polyurethane coating.
- Fire standard: EN 597-1, -2.



Olivia H^(1,2)

- Colour: grey.
- Zip on three sides.
- Two transportation handles on right-hand side.
- Composition: 55 % Polyester/45 % Polyurethane coating.
- Fire standard: EN 597-1, -2.



Marin^(1,2)

- Colour: blue
- Zip on one or three sides.
- Composition: 37 % Polyester/63 % Polyurethane coating.
- Fire standard: EN-597-1, -2.



Marin FR^(1,2)

- Colour: blue
- Zip on one short side.
- Composition: 55 % Polyester/45 % Polyurethane coating.
- Extra flame retardant, SS 876 00 10.



Nordic^(1,2)

- Colour: blue, grey.
- Zip on three sides.
- Composition: 44 % Polyamide/56 % Polyurethane coating.
- Welded seams.
- Handle on one long side.
- Fire standard: EN-597-1, -2.



Skandic^(1,2)

- Colour: grey, blue.
- Moisture-proof zip on one short side.
- Welded seams.
- Composition: 44 % Polyamide/56 % Polyurethane coating.
- Fire standard: EN-597-1, -2.



Steel^(1,2)

- Colour: blue
- Drip-proof zip on four sides.
- Welded seams.
- Composition: 44 % Polyamide/56 % Polyurethane coating.
- Fire standard: EN-597-1, -2.



Stone^(1,2)

- Colour: grey.
- Drip-proof zip on four sides.
- Welded seams.
- Composition: 61 % Polyamide/39 % Polyurethane coating.
- Fire standard: EN-597-1, -2, BS 6807 Crib 5.



BOB Extreme/Black^(1,2)

- Colour: blue, black.
- Two sections, zip on one short side.
- Handles on long and short sides.
- Suitable for Optimal BOB mattress.
- Composition: 40 % Polyamide/60 % Polyurethane.
- Fire standard: EN-597-1, -2.

Facts

- Moisture barrier, vapour permeable, vapour is removed, reducing risk of skin maceration.
- Soft, elastic.
- Simple care/handling.
- Wipe-clean, surface disinfection, wash at 95°C, can be tumble dried.
- Available in different sizes.
- Drip-proof zip.
- Complies with SS 876 00 1.

www.careofsweden.se

Bottom parts

Bottom Part Evac^(1,2)

- Colour: black.
- Zip, drip-proof, open-ended.
- Velcro strips for securing in emergency situation.
- Functional handles for movement/transport.
- Composition: 100 % Polyester Polyurethane.
- Can be combined with different top parts.
- Fire standard: EN-597-1, -2, BS 6807 Crib 5.



Bottom Part CuroCell[®]^(1,2)

- Colour: black.
- Composition: 100 % Polyester Polyurethane.
- Fire standard: EN-597-1, -2, BS 6807 Crib 5.



Bottom Part CuroCell[®] Extreme^(1,2)

- Colour: black.
- Composition: 40 % Polyamide/60 % Polyurethane coating.
- Fire standard: EN-597-1, -2, BS 6807 Crib 5.



Bottom Part Duro^(1,2)

- Colour: blue
- Drip-proof zip on four sides.
- Welded seams.
- Composition: 44 % Polyamide/56 % Polyurethane coating.
- Fire standard: EN-597-1, -2.



Optimal 5zon Plus[®] Bottom Part^(1,2)

- Colour: black.
- Functional handles for movement of mattress, short and long sides.
- Suitable for Optimal 5zon Plus mattress.
- Composition: 100 % Polyester Polyurethane.
- Fire standard: EN-597-1, -2.



Other mattress covers (not hygiene)

Quilt

- Colour: grey.
- Cotton cover.
- Wash at 70°C.
- Composition: 37 % Cotton/34 % Polyamide/29 % Polyester.



Triåk

- Knitted cotton.
- Colour: white.
- Washable at 60 °C, can be tumble dried.
- Composition: 100 % Cotton.



References

(1) ISO 16603 "resistance to penetration by blood and body fluids", ISO 16604 "resistance to penetration by blood-borne pathogens", EN 14126, SS 876 00 19 "Health care textiles – Bacteria penetration – Wet".

(2) SS-EN ISO 15496:2004, DIN 53.122.



Read the user manual prior to using the product.

Care of Sweden, CuroCell and Optimal 5zon are registered trademarks of Care of Sweden AB. © Care of Sweden AB, 2017. All rights reserved. Last updated: 2017-02-07.