



Mattress Covers

Meets health care requirements
for hygiene, quality and safety

Hygiene covers (full cover or top part only)



Olivia^(1,2,3)

- Colour: grey.
- Zip on three sides.
- Composition: 55 % Polyester/45 % Polyurethane coating.
- Fire standard: EN 597-1, -2.



Olivia H^(1,2,3)

- Colour: grey.
- Zip on three sides.
- Two transportation handles on right-hand side.
- Composition: 55 % Polyester/45 % Polyurethane coating.
- Fire standard: EN 597-1, -2.



Marin^(1,2,3)

- Colour: blue
- Zip on one or three sides.
- Composition: 37 % Polyester/63 % Polyurethane coating.
- Fire standard: EN-597-1, -2.



Marin FR^(1,2,3)

- Colour: blue
- Zip on one short side.
- Composition: 55 % Polyester/45 % Polyurethane coating.
- Extra flame retardant, SS 876 00 10.



Nordic^(1,2,3)

- Colour: blue, grey.
- Zip on three sides.
- Composition: 44 % Polyamide/56 % Polyurethane coating.
- Welded seams.
- Handle on one long side.
- Fire standard: EN-597-1, -2.



Skandic^(1,2,3)

- Colour: grey, blue.
- Moisture-proof zip on one short side.
- Welded seams.
- Composition: 44 % Polyamide/56 % Polyurethane coating.
- Fire standard: EN-597-1, -2.



Steel^(1,2,3)

- Colour: blue
- Drip-proof zip on four sides.
- Welded seams.
- Composition: 44 % Polyamide/56 % Polyurethane coating.
- Fire standard: EN-597-1, -2.



Stone^(1,2,3)

- Colour: grey.
- Drip-proof zip on four sides.
- Welded seams.
- Composition: 61 % Polyamide/39 % Polyurethane coating.
- Fire standard: EN-597-1, -2, BS 6807 Crib 5.



BOB Extreme/Black^(1,2,3)

- Colour: blue, black.
- Two sections, zip on one short side.
- Handles on long and short sides.
- Suitable for Optimal BOB mattress.
- Composition: 40 % Polyamide/60 % Polyurethane.
- Fire standard: EN-597-1, -2.

Facts

- Moisture barrier, vapour permeable, vapour is removed, reducing risk of skin maceration.
- Elastic.
- Simple care/handling.
- Wipe-clean, surface disinfection, wash at 95°C, can be tumble dried.
- Available in different sizes.
- Drip-proof zip.
- Complies with SS 876 00 1.

Bottom parts

Bottom Part Evac^(1,2,3)

- Colour: black.
- Zip, drip-proof, open-ended.
- Velcro strips for securing in emergency situation.
- Functional handles for movement/transport.
- Composition: 100 % Polyester Polyurethane.
- Can be combined with different top parts.
- Fire standard: EN-597-1, -2, BS 6807 Crib 5.



Bottom Part CuroCell^{®(1,2,3)}

- Colour: black.
- Composition: 100 % Polyester Polyurethane.
- Fire standard: EN-597-1, -2, BS 6807 Crib 5.



Bottom Part CuroCell^{® Extreme(1,2,3)}

- Colour: black.
- Composition: 40 % Polyamide/60 % Polyurethane coating.
- Fire standard: EN-597-1, -2, BS 6807 Crib 5.



Bottom Part Duro^(1,2,3)

- Colour: blue
- Drip-proof zip on four sides.
- Welded seams.
- Composition: 44 % Polyamide/56 % Polyurethane coating.
- Fire standard: EN-597-1, -2.



Optimal 5zon Plus^{® Bottom Part(1,2,3)}

- Colour: black.
- Functional handles for movement of mattress, short and long sides.
- Suitable for Optimal 5zon Plus mattress.
- Composition: 100 % Polyester Polyurethane.
- Fire standard: EN-597-1, -2.



Other mattress covers (not hygiene)

Quilt

- Colour: grey.
- Cotton cover.
- Wash at 70°C.
- Composition: 37 % Cotton/34 % Polyamide/29 % Polyester.



Triåk

- Knitted cotton.
- Colour: white.
- Washable at 60 °C, can be tumble dried.
- Composition: 100 % Cotton.



References

(1) CDC (2003), Guidelines for Environmental Infection Control in Health-Care facilities, updated 2019.

(2) SS-EN ISO 15496:2004, DIN 53.122.

(3) SS 876 00 20:2017, Healthcare textiles - Mattresses - Specifications and requirements



Read the instruction for use prior to using the product.

Care of Sweden, CuroCell and Optimal 5zon are registered trademarks of Care of Sweden AB. © Care of Sweden AB, 2017. All rights reserved.
Last updated: 2021-06-08.