

# PRIMA 6



A functional elevating pad

# PRIMA 6

PRIMA 6 is an easy and functional basic mattress with good durability.

It is an elevating pad that is placed underneath a mattress overlay, for example CuroCell® Nova (air mattress) or Optimal 5zon® BM (preventative).

Care of Sweden's mattresses are supplied with a four way stretch hygienic outer cover which conforms to relevant standards for hygiene<sup>(1)</sup>, quality and safety. All covers are manufactured from modern, recyclable materials and are vapour permeable<sup>(2)</sup>, which reduce the risk of skin maceration.

## Facts

Core: HR foam 35 kg/m<sup>3</sup>

Height: 6 cm

Fire: EN 597-1, -2

Warranty: 2 years

Cleaning of cover: Wipe with cleaning agent and/or disinfectants.

Machine wash max 95 °C, tumble drying.



## Features & Characteristics

- Easy to clean
- Hygienic cover
- Reversible
- Unique labelling

## Available covers

- Marin with zip on one side
- Marin with zip on three sides.

## Sizes

80/85/90/105/120x200/210 cm.

If you need any other sizes please contact your distributor.

## References

(1) ISO 16603 "resistance of penetration by blood and body fluids", ISO 16604 "resistance of penetration by blood-borne pathogens", EN 14126, SS 876 00 19 "Health care textiles – Bacteria penetration – Wet".

(2) SS-EN ISO 15496:2004, DIN 53.122.



Always read the instruction for use prior to usage of the product.

Care of Sweden, CuroCell and Optimal 5zon are trademarks of Care of Sweden AB.  
© Care of Sweden AB, 2014. All rights reserved. Last update: 2014-03-17.