

PRIMA 6

A functional elevating pad



PRIMA 6 with cover Marin

PRIMA 6

PRIMA 6 is an easy and functional basic mattress with good durability.

It is an elevating pad that is placed underneath a mattress overlay, for example CuroCell® Nova (air mattress) or Optimal 5zon® BM (preventative).

Care of Sweden's mattresses are supplied with a four way stretch hygienic outer cover which conforms to relevant standards for hygiene, quality and safety⁽¹⁾⁽³⁾.

All covers are manufactured from modern, recyclable materials and are vapour permeable⁽²⁾, which reduce the risk of skin maceration.



AVAILABLE COVERS

Marin

With zipper on one side.

Marin

With zipper on three sides.

REFERENCES

(1) CDC (2003), Guidelines for Environmental Infection Control in Health-Care facilities, updated 2019

(2) SS-EN ISO 15496:2004, DIN 53122-1.

(3) SS 876 00 20:2017, Healthcare textiles - Mattresses - Specifications and requirements



Always read the instruction for use prior to usage of the product.

Care of Sweden, CuroCell and Optimal 5zon are trademarks of Care of Sweden AB. © Care of Sweden AB, 2021. All rights reserved. Last update: 2021-05-28.

FEATURES & CHARACTERISTICS	FACTS
<ul style="list-style-type: none">• Easy to clean• Hygienic cover• Reversible	<ul style="list-style-type: none">• Core: HR foam 35 kg/m³• Height: 6 cm• Fire: EN 597-1, -2• Warranty: 2 years• Cleaning of cover: Wipe with cleaning agent and/or disinfectants• Machine wash max 95 °C, tumble drying
SIZES	
80/85/90/105/120x200/210 cm. If you need any other sizes please contact your distributor.	