


PRIMA 12

Functional mattress replacement that offers
good comfort & durability



PRIMA 12 with cover Marin/Evac

PRIMA 12

PRIMA 12 has a profiled mattress core that provides user good comfort and pressure relief.

The mattress is delivered with a detachable cover Marin that has a zipper on one short side. As an option, you can choose the Marin/Evac cover which is detachable and has a drip shield zip and feature handles for evacuation.

PRIMA 12 is placed directly on the bed.

Care of Sweden's mattresses are supplied with a four way stretch hygienic outer cover which conforms to relevant standards for hygiene, quality and safety⁽¹⁾⁽³⁾.

All covers are manufactured from modern, recyclable materials and are vapour permeable⁽²⁾, which reduce the risk of skin maceration.



AVAILABLE COVERS

Marin

Drip shield zip on one side.

Marin

Drip shield zip on three sides.

Marin

Drip shield zip on three sides and handles on one side.

Marin/Evac

Divisible, drip shield zip and handles for evacuation on one side.

REFERENCES

(1) CDC (2003), Guidelines for Environmental Infection Control in Health-Care facilities, updated 2019

(2) SS-EN ISO 15496:2004, DIN 53.122.

(3) SS 876 00 20:2017, Healthcare textiles - Mattresses - Specifications and requirements.



Always read the instruction for use prior to usage of the product.

Care of Sweden is a trademark of Care of Sweden AB.
© Care of Sweden AB, 2021. All rights reserved.
Last update: 2021-05-28.

FEATURES & CHARACTERISTICS	FACTS
<ul style="list-style-type: none">• Easy to clean• Hygienic cover• Reversible	<ul style="list-style-type: none">• Core: HR foam 35 kg/m³• Height: 13 cm• Recommended user weight: up to 130 kg• Fire: EN 597-1, -2• Warranty: 3 years• Cleaning of cover: Wipe with cleaning agent and/or disinfectants• Machine wash max 95 °C, tumble drying.
SIZES	
80/85/90/105/120x200 cm. If you need any other sizes please contact your distributor.	