

PRIMA 12



The picture shows a mattress with Marin cover.

Functional mattress replacement
that offers good comfort & durability

PRIMA 12

PRIMA 12 has a profiled mattress core that provides user good comfort and pressure relief.

The mattress is delivered with a detachable cover Marin that has a zipper on one short side. As an option, you can choose

the Marin/Evac cover which is detachable and has a drip shield zip and feature handles for evacuation.

PRIMA 12 is placed directly on the bed.

Care of Sweden's mattresses are supplied with a four way stretch hygienic outer cover which conforms to relevant standards for hygiene⁽¹⁾, quality and safety. All covers are manufactured from modern, recyclable materials and are vapour permeable⁽²⁾, which reduce the risk of skin maceration.

Facts

Core: HR foam 35 kg/m³

Height: 13 cm

Recommended user weight: 0-130 kg

Fire: EN 597-1, -2

Warranty: 3 years

Cleaning of cover: Wipe with cleaning agent and/or disinfectants.

Machine wash max 95 °C, tumble drying.



Features & Characteristics

- Easy to clean
- Hygienic cover
- Reversible
- Unique labelling

Available covers

Marin

Drip shield zip on one side.

Marin

Drip shield zip on three sides.

Marin

Drip shield zip on three sides and handles on one side.

Marin/Evac

Divisible, drip shield zip and handles for evacuation on one side.

In need of non hygienic cover, please contact your distributor.

Sizes

80/85/90/105/120x200 cm.

If you need any other sizes please contact your distributor.

References

(1) ISO 16603 "resistance of penetration by blood and body fluids", ISO 16604 "resistance of penetration by blood-borne pathogens", EN 14126, SS 876 00 19 "Health care textiles – Bacteria penetration – Wet".

(2) SS-EN ISO 15496:2004, DIN 53.122.



Always read the instruction for use prior to usage of the product.

Care of Sweden is a trademark of Care of Sweden AB.

© Care of Sweden AB, 2017. All rights reserved. Last update: 2017-05-03.



Mattress selector

- Guidelines when choosing a product
See how to choose the right mattress at:
www.careofsweden.com/mattress-selector

 **care of sweden**
Supporting Life

Care of Sweden AB, Box 146, SE-514 23 Tranemo, Sweden
+46 (0)771 106 600 | info@careofsweden.se | www.careofsweden.com