


# PRIMA 12



The picture shows a mattress with Marin/Evac cover.

Functional mattress replacement  
that offers good comfort & durability

# PRIMA 12

PRIMA 12 has a profiled mattress core that provides user good comfort and pressure relief.

The mattress is delivered with a detachable cover Marin that has a zipper on one short side. As an option, you can choose

the Marin/Evac cover which is detachable and has a drip shield zip and feature handles for evacuation.

PRIMA 12 is placed directly on the bed.

Care of Sweden's mattresses are supplied with a four way stretch hygienic outer cover which conforms to relevant standards for hygiene<sup>(1)</sup>, quality and safety. All covers are manufactured from modern, recyclable materials and are vapour permeable<sup>(2)</sup>, which reduce the risk of skin maceration.

## Facts

Core: HR foam 35 kg/m<sup>3</sup>

Height: 13 cm

Recommended user weight: 0-130 kg

Fire: EN 597-1, -2

Warranty: 3 years

Cleaning of cover: Wipe with cleaning agent and/or disinfectants.

Machine wash max 95 °C, tumble drying.



## Features & Characteristics

- Easy to clean
- Hygienic cover
- Reversible
- Unique labelling

## Available covers

### Marin

Drip shield zip on one side.

### Marin

Drip shield zip on three sides.

### Marin

Drip shield zip on three sides and handles on one side.

### Marin/Evac

Divisible, drip shield zip and handles for evacuation on one side.

In need of non hygienic cover, please contact your distributor.

## Sizes

80/85/90/105/120x200 cm.

If you need any other sizes please contact your distributor.

## References

(1) ISO 16603 "resistance of penetration by blood and body fluids", ISO 16604 "resistance of penetration by blood-borne pathogens", EN 14126, SS 876 00 19 "Health care textiles – Bacteria penetration – Wet".

(2) SS-EN ISO 15496:2004, DIN 53.122.



Always read the instruction for use prior to usage of the product.

Care of Sweden is a trademark of Care of Sweden AB.

© Care of Sweden AB, 2017. All rights reserved. Last update: 2017-05-03.



### Mattress selector

- Guidelines when choosing a product  
See how to choose the right mattress at:  
[www.careofsweden.com/mattress-selector](http://www.careofsweden.com/mattress-selector)

 **care of sweden**  
Supporting Life

Care of Sweden AB, Box 146, SE-514 23 Tranemo, Sweden  
+46 (0)771 106 600 | [info@careofsweden.se](mailto:info@careofsweden.se) | [www.careofsweden.com](http://www.careofsweden.com)