

PRIMA BM



Mattress overlay
that offers good comfort & durability

PRIMA BM

Mattress overlay PRIMA BM has a mattress core in one piece that provides the patient good comfort and pressure relief. It is delivered with hygienic cover Marin.

Care of Sweden's mattresses are supplied with a four way stretch hygienic outer cover which conforms to relevant standards for hygiene⁽¹⁾, quality and safety. All covers are manufactured from modern, recyclable materials and are vapour permeable⁽²⁾, which reduce the risk of skin maceration.

Facts

Core: HR foam 35 kg/m³

Height: 6 cm

Recommended user weight: 0-130 kg*

Fire: EN 597-1, -2

Warranty: 2 years

Cleaning of cover: Wipe with cleaning agent and/or disinfectants.

Machine wash max 95 °C, tumble drying.

*Please note that the underlying mattress affects the overall function including the recommended user weight.



Features & Characteristics

- Easy to clean
- Hygienic cover
- Reversible
- Unique labelling

Available covers

- Marin with zip on one side
- Marin with zip on three sides.

In need of non hygienic cover, please contact your distributor.

Sizes

80/85/90/105/120x200 cm.

If you need any other sizes please contact your distributor.

References

(1) ISO 16603 "resistance of penetration by blood and body fluids", ISO 16604 "resistance of penetration by blood-borne pathogens", EN 14126, SS 876 00 19 "Health care textiles – Bacteria penetration – Wet".

(2) SS-EN ISO 15496:2004, DIN 53.122.



Always read the instruction for use prior to usage of the product.

Care of Sweden is a trademark of Care of Sweden AB.

© Care of Sweden AB, 2015. All rights reserved. Last update: 2015-02-12.



Mattress selector

- Guidelines when choosing a product
See how to choose the right mattress at:
www.careofsweden.com/mattress-selector

 care of sweden
Supporting Life

Care of Sweden AB, Box 146, SE-514 23 Tranemo, Sweden
+46 (0)771 106 600 | info@careofsweden.se | www.careofsweden.com