

Optimal M3



The picture shows a mattress with Marin cover.

High quality mattress replacement

Optimal M3

Optimal M3 is our best-selling and much-loved standard mattress, used in all types of healthcare environments.

It is a split core mattress with high density and is delivered with hygienic cover Marin and Marin/Evac. The Marin cover is detachable with a zip fastener on one of its short sides. The

Marin/Evac cover is equipped with a divisible, drip shield zip and feature handles for evacuation.

Optimal M3 is a high quality and economical solution for everyday use.

Care of Sweden's mattresses are supplied with a four way stretch hygienic outer cover which conforms to relevant standards for hygiene⁽¹⁾, quality and safety. All covers are manufactured from modern, recyclable materials and are vapour permeable⁽²⁾, which reduce the risk of skin maceration.

Facts

Core: HR foam 50 kg/m³

Height: 12 cm

Recommended user weight: 0-130 kg

Fire: EN 597-1, -2

Warranty: 3 years

Cleaning of cover: Wipe with cleaning agent and/or disinfectants.

Machine wash max 95 °C, tumble drying.



Features & Characteristics

Pressure

Highly elastic, high-quality material with effective pressure equalisation that provides good comfort.

Shear forces

The mattress is specially designed to minimise the effects of shear forces.

Moisture

The hygienic cover is vapour permeable⁽²⁾, which reduces the risk of skin maceration.

Hygiene

The cover is easy to clean and has a very well documented hygienic standard⁽¹⁾.

Available covers

Marin

Drip shield zip on one side.

Marin

Drip shield zip on three sides.

Marin

Drip shield zip on three sides and handles on one side.

Marin/Evac

Divisible, drip shield zip and handles for evacuation on one side.

In need of non hygienic cover, please contact your distributor.

Sizes

80/85/90/105/120x200 cm and 80/85/90x205/210 cm.

If you need any other sizes please contact your distributor.

References

(1) ISO 16603 "resistance of penetration by blood and body fluids", ISO 16604 "resistance of penetration by blood-borne pathogens", EN 14126, SS 876 00 19 "Health care textiles – Bacteria penetration – Wet".

(2) SS-EN ISO 15496:2004, DIN 53.122.



Always read the instruction for use prior to usage of the product.

Care of Sweden is a trademark of Care of Sweden AB.

© Care of Sweden AB, 2017. All rights reserved. Last update: 2017-05-03.



Mattress selector

- Guidelines when choosing a product
See how to choose the right mattress at:
www.careofsweden.com/mattress-selector

 **care of sweden**
Supporting Life

Care of Sweden AB, Box 146, SE-514 23 Tranemo, Sweden
+46 (0)771 106 600 | info@careofsweden.se | www.careofsweden.com