

# Optimal 5zon<sup>®</sup>



The picture shows a mattress with Stone/Evac cover.

The complete mattress  
for pressure ulcer prevention

# Optimal 5zon<sup>®</sup>

The Optimal 5zon<sup>®</sup> is ideal for all healthcare environments and is used as an aid in the prevention and treatment of pressure ulcers up to and including category 2<sup>(1)</sup>, and for pain therapy. It is a soft and flexible mattress replacement that is placed directly on the bed and replaces both the standard mattress and mattress pad.

It has very good pressure redistribution properties<sup>(2)</sup> and is also equipped with the unique DuoCore<sup>®</sup> sit function. It has been designed to provide support when seated, while also offering good comfort in lying position.

Optimal 5zon<sup>®</sup> comes with several different options of covers.

Care of Sweden's mattresses are supplied with a four way stretch hygienic outer cover which conforms to relevant standards for hygiene<sup>(3)</sup>, quality and safety. All covers are manufactured from modern, recyclable materials and are vapour permeable<sup>(4)</sup>, which reduce the risk of skin maceration.

## Facts

Core: HR-foam 65+50 kg/m<sup>3</sup>

Height: 14 cm

Function: Prevention and treatment up to and including category 2<sup>(1)</sup>

Recommended user weight: 0-180 kg

Fire: EN 597-1, -2

Warranty: 3 years

Cleaning of cover: Wipe with cleaning agent and/or disinfectants.

Machine wash max 95 °C, tumble drying.



## Features & Characteristics

### Pressure

Highly elastic, high-quality material and a unique zone and module system gives effective pressure equalisation.

### Shear forces

The mattress is specially designed to minimise the effects of shear forces.

### Sit function

Our patented<sup>(5)</sup> sit function DuoCore<sup>®</sup> is designed to give support when sitting and also to provide good comfort when lying.

### Moisture

The hygienic cover is vapour permeable<sup>(4)</sup>, which reduces the risk of skin maceration.

### Hygiene

The cover is easy to clean and has a very well documented hygienic standard<sup>(3)</sup>.

## Available covers

### Olivia/Evac Cover

2-piece cover with handles for evacuation. Low friction top cover is also available.

### Cover Olivia

Drip shield zip on three sides.

### Cover Olivia H

Same as Olivia cover but with handles for transport and moving.

### Cover Nordic

Welded seams, drip shield zip on three sides and handles on one side.

### Steel/Evac Cover

2-piece cover with handles for evacuation. Welded seams and drip shield zip on all sides.

### Cover Stone

Crib 5 cover with welded seams and a protected zipper on three sides.

### Stone/Evac Cover

2-piece Crib 5 cover with handles for evacuation. Welded seams and drip shield zip on all sides.

For need of non hygienic cover, please contact your distributor.

## Sizes

80/85/90/105/120x200 cm and 80/85/90x210 cm.

If you need any other sizes please contact your distributor.

## References

(1) National Pressure Ulcer Advisory Panel, European Pressure Ulcer Advisory Panel and Pan Pacific Pressure Injury Alliance. Prevention and Treatment of Pressure Ulcers: Quick Reference Guide. Emily Haesler (Ed.). Cambridge Media: Osborne Park, Western Australia; 2014.

(2) Examination of the shear force and the pressure relief capability, Berlincert 2009.

(3) ISO 16603 "resistance of penetration by blood and body fluids", ISO 16604 "resistance of penetration by blood-borne pathogens", EN 14126, SS 876 00 19 "Health care textiles – Bacteria penetration – Wet".

(4) SS-EN ISO 15496:2004, DIN 53.122.

(5) Patent no. 008522741.



Always read the instruction for use prior to usage of the product.

Care of Sweden, Optimal 5zon and DuoCore are trademarks of Care of Sweden AB. © Care of Sweden AB, 2017. All rights reserved. Last update: 2017-09-26.



## Mattress selector

- Guidelines when choosing a product  
See how to choose the right mattress at:  
[www.careofsweden.com/mattress-selector](http://www.careofsweden.com/mattress-selector)

 care of sweden  
Supporting Life

Care of Sweden AB, Box 146, SE-514 23 Tranemo, Sweden  
+46 (0)771 106 600 | [info@careofsweden.se](mailto:info@careofsweden.se) | [www.careofsweden.com](http://www.careofsweden.com)