

Optimal 5zon[®] BM

Preventative mattress overlay
with good comfort



Optimal 5zon BM[®] with cover Olivia

Optimal 5zon[®] BM

The Optimal 5zon[®] BM is a flexible mattress overlay that is intended for community and home care, as an aid in helping to prevent and treat pressure ulcers up to category 2⁽¹⁾.

Its design provides good pressure redistribution and, thanks to the High Resilience Foam core, the Optimal 5zon[®] BM has a long life and provides a good level of comfort.

The cover is also equipped with a drip-shield zip on three sides and elasticated fastening straps.

Care of Sweden's mattresses are supplied with a four way stretch hygienic outer cover which conforms to relevant standards for hygiene, quality and safety⁽²⁾⁽⁴⁾. All covers are manufactured from modern, recyclable materials and are vapour permeable⁽³⁾, which contributes to reduce the risk of skin maceration.



SIZES

80/85/90x200/205/210 cm and 105/120x200/210 cm.

If you need any other sizes please contact your distributor.

REFERENCES

(1) European Pressure Ulcer Advisory Panel, National Pressure Injury Advisory Panel and Pan Pacific Pressure Injury Alliance. Prevention and Treatment of Pressure Ulcers/Injuries: Clinical Practice Guideline. The International Guideline. Emily Haesler (Ed). EPUAP/NPIAP/PPPIA: 2019

(2) CDC (2003), Guidelines for Environmental Infection Control in Health-Care facilities, updated 2019.

(3) SS-EN ISO 15496:2004, DIN 53122-1.

(4) SS 876 00 20:2017, Healthcare textiles - Mattresses - Specifications and requirements



Always read the instructions for use prior to usage of the product.

Care of Sweden, Optimal 5zon and DuoCore are trademarks of Care of Sweden AB. © Care of Sweden AB, 2021. All rights reserved. Last update: 2021-06-03.

FEATURES & CHARACTERISTICS	FACTS
<p>Pressure Elastic material and a zone and module system gives good pressure equalisation.</p> <p>Shear forces The mattress is specially designed to reduce the effects of shear forces.</p> <p>Moisture The hygienic cover is vapour permeable⁽³⁾, which contributes to reduce the risk of skin maceration.</p> <p>Hygiene The cover is easy to clean and has a good hygienic standard⁽²⁾⁽⁴⁾.</p>	<ul style="list-style-type: none"> • Core: HR-foam 65 kg/m³ • Height: 6 cm • Covers: Olivia • Function: as an aid for prevention and treatment up to and including category 2⁽¹⁾ • Recommended user weight: up to 150 kg* • Fire: EN 597-1, -2 • Warranty: 2 years • Cleaning of cover: wipe with cleaning agent and/or disinfectants. • Machine wash max 95 °C, tumble drying. <p>*Please note that the underlying mattress affects the overall function including the recommended user weight.</p>