

Optimal Care

Pressure redistribution mattress with good comfort



Optimal Care with cover Olivia/Evac

FACTS	SIZES
<ul style="list-style-type: none"> • Core: HR foam 50 kg/m³ • Height: 14 cm • Function: an aid for prevention and treatment up to and including category 2⁽¹⁾ • Recommended user weight: up to 160 kg • Fire: EN 597-1, -2 • Warranty: 3 years • Cleaning of cover: wipe with cleaning agent and/or disinfectants. • Machine wash max 95 °C, tumble drying. 	<p>80/85/90x200/205/210 cm and 105/120x200 cm. If you need any other sizes please contact your distributor.</p>

Optimal Care

A modern mattress should have good pressure ulcer prevention properties, provide good comfort and be easy to handle. It should also represent good health economics. Optimal Care has been made to meet the requirements for a modern mattress used in health care environments.

Optimal Care is a pressure redistributing mattress divided into zones adapted to the body. It also has reinforced sides to give support when getting into and out of bed. Optimal Care is a mattress designed to meet stringent requirements in terms of quality, comfort and support.

The label contains the necessary information to ensure easy use of this mattress. Optimal Care comes with several different options of covers.

Care of Sweden's mattresses are supplied with a four way stretch hygienic outer cover which conforms to relevant standards for hygiene, quality and safety⁽²⁾⁽⁴⁾. All covers are manufactured from modern, recyclable materials and are vapour permeable⁽³⁾, which contributes to reduce the risk of skin maceration.

REFERENCES

(1) European Pressure Ulcer Advisory Panel, National Pressure Injury Advisory Panel and Pan Pacific Pressure Injury Alliance. Prevention and Treatment of Pressure Ulcers/Injuries: Clinical Practice Guideline. The International Guideline. Emily Haesler (Ed). EPUAP/NPIAP/PPPIA: 2019

(2) CDC (2003), Guidelines for Environmental Infection Control in Health-Care facilities, updated 2019.

(3) SS-EN ISO 15496:2004, DIN 53.122.

(4) SS 876 00 20:2017, Healthcare textiles - Mattresses - Specifications and requirements



FEATURES & CHARACTERISTICS	AVAILABLE COVERS
<p>Pressure Elastic material with a module system gives pressure equalisation.</p> <p>Shear forces The mattress is specially designed to reduce the effects of shear forces.</p> <p>Moisture The hygienic cover is vapour permeable⁽³⁾, which contributes to reduce the risk of skin maceration.</p> <p>Hygiene The cover is easy to clean and has a good hygienic standard⁽²⁾⁽⁴⁾.</p>	<p>Olivia Drip shield zip on three sides.</p> <p>Olivia H Same as Olivia cover but with handles for transport and moving.</p> <p>Evac 2-piece cover with handles for evacuation. Low friction top cover is also available.</p> <p>Nordic Welded seams, drip shield zip on three sides and handles on one side.</p> <p> Always read the instructions for use prior to usage of the product.</p> <p>Care of Sweden is a trademark of Care of Sweden AB. © Care of Sweden AB, 2021. All rights reserved. Last update: 2021-06-03.</p>