

Optimal Cura

Preventative mattress replacement with
good pressure equalization



Optimal Cura with cover Olivia

Optimal Cura

Optimal Cura is used as an aid in the prevention and treatment of pressure ulcers up to and including category 2⁽¹⁾. It is a mattress replacement, used in all types of healthcare environments.

Optimal Cura is a two layer mattress where the top layer is made from visco-elastic foam (so called memory foam) and the bottom layer is made from HR foam. The design of the mattress and the materials gives good pressure equalisation. The two layers allows the body to sink into the material, and thereby lower the pressure against the body.

Care of Sweden's mattresses are supplied with a four way stretch hygienic outer cover which conforms to relevant standards for hygiene, quality and safety⁽²⁾⁽⁴⁾. All covers are manufactured from modern, recyclable materials and are vapour permeable⁽³⁾, which reduce the risk of skin maceration.



AVAILABLE COVERS

Cover Olivia

Drip shield zip on three sides.

Cover Olivia H

Same as Olivia cover but with handles for transport and moving.

Cover Olivia/Evac

Divisible, drip shield zip and handles for evacuation on one side.

REFERENCES

(1) European Pressure Ulcer Advisory Panel, National Pressure Injury Advisory Panel and Pan Pacific Pressure Injury Alliance. Prevention and Treatment of Pressure Ulcers/Injuries: Clinical Practice Guideline. The International Guideline. Emily Haesler (Ed). EPUAP/NPIAP/PPPIA: 2019

(2) CDC (2003), Guidelines for Environmental Infection Control in Health-Care facilities, updated 2019.

(3) SS-EN ISO 15496:2004, DIN 53122-1.

(4) SS 876 00 20:2017, Healthcare textiles - Mattresses - Specifications and requirements



Always read the instruction for use prior to usage of the product.

Care of Sweden are trademarks of Care of Sweden AB. ©

Care of Sweden AB, 2021. All rights reserved. Last update: 2021-07-06.

FEATURES & CHARACTERISTICS	FACTS
<p>Pressure Elastic material with good pressure equalisation that provides good comfort.</p> <p>Shear forces The mattress is designed to reduce the effects of shear forces.</p> <p>Moisture The hygienic cover is vapour permeable⁽³⁾, which contributes to reduce the risk of skin maceration.</p> <p>Hygiene The cover is easy to clean and has a good documented hygienic standard⁽²⁾⁽⁴⁾.</p>	<ul style="list-style-type: none"> • Core: HR-foam + visco-elastic foam 50 kg/m³ • Height: 13 cm • Function: Prevention and treatment up to and including category 2⁽¹⁾ • Recommended user weight: up to 130 kg • Fire: EN 597-1, -2 • Warranty: 3 years • Cleaning of cover: Wipe with cleaning agent and/or disinfectants. • Machine wash max 95 °C, tumble drying.
	<p>SIZE</p> <p>80/85/90/105/120x200 cm and 80/85/90x205/210 cm. If you need any other sizes please contact your distributor.</p>