

Optimal Solett

Preventative mattress replacement for
lightweight and sensitive patients



Optimal Solett with Olivia/Evac cover

Optimal Solett

Optimal Solett is used as an aid in the prevention and treatment of pressure ulcers up to category 1⁽¹⁾. It has been specially designed for lightweight and sensitive users. This mattress replacement provides good pressure redistribution.

Optimal Solett is delivered with the Olivia/Evac cover. It is equipped with a divisible, drip shield zip and feature handles for evacuation.

Care of Sweden's mattresses are supplied with a four way stretch hygienic outer cover which conforms to relevant standards for hygiene, quality and safety⁽²⁾⁽⁴⁾. All covers are manufactured from modern, recyclable materials and are vapour permeable⁽³⁾, which reduce the risk of skin maceration.

SIZES

80/85/90/105/120x200 cm and 80/85/90x210 cm. If you need any other sizes please contact your distributor.



AVAILABLE COVERS

Olivia/Evac Cover

Divisible, drip shield zip and handles for evacuation on one side.

Cover Olivia

Drip shield zip on three sides.

Cover Nordic

Welded seams, drip shield zip on three sides and handles on one side.

Steel/Evac Cover

2-piece cover with handles for evacuation. Welded seams and a protected zipper on all sides.

For need of non hygienic cover, please contact your distributor.


REFERENCES

(1) European Pressure Ulcer Advisory Panel, National Pressure Injury Advisory Panel and Pan Pacific Pressure Injury Alliance. Prevention and Treatment of Pressure Ulcers/Injuries: Clinical Practice Guideline. The International Guideline. Emily Haesler (Ed). EPUAP/NPIAP/PPPIA: 2019

(2) CDC (2003), Guidelines for Environmental Infection Control in Health-Care facilities, updated 2019

(3) SS-EN ISO 15496:2004, DIN 53122-1

(4) SS 876 00 20:2017, Healthcare textiles - Mattresses - Specifications and requirements

FEATURES & CHARACTERISTICS	FACTS
<p>Pressure Elastic material with pressure equalisation that provides good comfort.</p> <p>Shear forces The mattress is specially designed to minimise the effects of shear forces.</p> <p>Moisture The hygienic cover is vapour permeable⁽³⁾, which reduces the risk of skin maceration.</p> <p>Hygiene The cover is easy to clean and has a good documented hygienic standard⁽²⁾⁽⁴⁾.</p>	<ul style="list-style-type: none"> • Core: HR-foam 35 kg/m³ + visco-elastic foam 50 kg/m³ • Height: 12 cm • Function: prevention and treatment up to and including category 1⁽¹⁾ • Recommended user weight: up to 60 kg • Fire: EN 597-1, -2 • Warranty: 2 years • Cleaning of cover: wipe with cleaning agent and/or disinfectants • Machine wash max 95 °C, tumble drying <p> Always read the instruction for use prior to usage of the product.</p> <p>Care of Sweden is a trademark of Care of Sweden AB. © Care of Sweden AB, 2021. All rights reserved. Last update: 2021-05-28.</p>