

Optimal Solett

Preventative mattress replacement for
lightweight and sensitive patients



Optimal Solett with Olivia/Evac cover

Optimal Solett

Optimal Solett is used as an aid in the prevention and treatment of pressure ulcers up to category 1⁽¹⁾. It has been specially designed for lightweight and sensitive patients. This mattress replacement provides very good pressure redistribution.

Optimal Solett is delivered with the Olivia/Evac cover. It is equipped with a divisible, drip shield zip and feature handles for evacuation.

Care of Sweden's mattresses are supplied with a four way stretch hygienic outer cover which conforms to relevant standards for hygiene⁽²⁾, quality and safety. All covers are manufactured from modern, recyclable materials and are vapour permeable⁽³⁾, which reduce the risk of skin maceration.

SIZES

80/85/90/105/120x200 cm and 80/85/90x210 cm.
If you need any other sizes please contact your distributor



AVAILABLE COVERS

Olivia/Evac Cover

Divisible, drip shield zip and handles for evacuation on one side.

Cover Olivia

Drip shield zip on three sides.

Cover Nordic

Welded seams, drip shield zip on three sides and handles on one side.

Steel/Evac Cover

2-piece cover with handles for evacuation. Welded seams and a protected zipper on all sides.

For need of non hygienic cover, please contact your distributor.


REFERENCES

- (1) Haesler, E. (Ed). National Pressure Ulcer Advisory Panel, European Pressure Ulcer Advisory Panel and Pan Pacific Pressure Injury Alliance. (2014). Prevention and Treatment of Pressure Ulcers: Clinical Guideline. Cambridge Media: Perth, Australia.
 (2) ISO 16603 "resistance of penetration by blood and body fluids", ISO 16604 "resistance of penetration by blood-borne pathogens", EN 14126.
 (3) SS-EN ISO 15496:2004, DIN 53122-1



Always read the instruction for use prior to usage of the product.

Care of Sweden is a trademark of Care of Sweden AB.
 © Care of Sweden AB, 2019. All rights reserved.
 Last update: 2019-08-19.

FEATURES & CHARACTERISTICS	FACTS
<p>Pressure Highly elastic, high-quality material with effective pressure equalisation that provides good comfort.</p> <p>Shear forces The mattress is specially designed to minimise the effects of shear forces.</p> <p>Moisture The hygienic cover is vapour permeable⁽³⁾, which reduces the risk of skin maceration.</p> <p>Hygiene The cover is easy to clean and has a very well documented hygienic standard⁽²⁾.</p>	<ul style="list-style-type: none"> • Core: HR-foam 35 kg/m³ + visco-elastic foam 50 kg/m³ • Height: 12 cm • Function: Prevention and treatment up to and including category 1⁽¹⁾ • Recommended user weight: 0-60 kg • Fire: EN 597-1, -2 • Warranty: 2 years • Cleaning of cover: Wipe with cleaning agent and/or disinfectants • Machine wash max 95 °C, tumble drying
	 <p>MATTRESS SELECTOR – Guidelines when choosing a product See how to choose the right mattress at: www.careofsweden.com/mattress-selector</p>