

CuroCell S.A.M.® CF16



Mattress replacement system
with good comfort

CuroCell S.A.M.[®] CF16

CuroCell S.A.M.[®] CF16 is used as an aid in the prevention and treatment of pressure ulcers up to category 2⁽¹⁾. It is a static air mattress which provides good comfort and pressure⁽²⁾ relief. It is a mattress replacement system with an integrated foam core.

The mattress is made of soft and durable materials. The included support pump is easy to use and quiet. The mattress comes with CPR (rapid deflation of air) and has a built-in heel function.

Thanks to the extra soft inner cover, the already low pressure⁽³⁾ is further reduced.

Care of Sweden's mattresses are supplied with a four way stretch hygienic outer cover which conforms to relevant standards for hygiene, quality and safety.⁽⁴⁾⁽⁶⁾ All covers are manufactured from modern, recyclable materials and are vapour permeable⁽⁵⁾, which reduce the risk of skin maceration.

Facts

Type: Static air mattress
Height: 16 cm
Function: Prevention and treatment up to and including category 2⁽¹⁾
Recommended user weight: 0-200 kg
Fire: EN 597-1, -2
Power supply: 230 V/50 Hz
Warranty: 2 years
Cleaning of cover: Wipe with cleaning agent and/or disinfectants.
Machine wash max 95 °C, tumble drying.



Features & Characteristics

- CPR function – rapid deflation
- Easy to handle and maintain
- Extra soft inner cover
- Heel function
- Silent and energy-efficient pump
- Soft, elastic, vapour permeable⁽⁵⁾ hygienic cover
- Good cleaning properties

Sizes

80/85/90/105/120x200/210 cm.
If you need any other sizes please contact your distributor.

Optional

- Durable bag for transportation and storage
- Transport function

References

- (1) European Pressure Ulcer Advisory Panel, National Pressure Injury Advisory Panel and Pan Pacific Pressure Injury Alliance. Prevention and Treatment of Pressure Ulcers/Injuries: Clinical Practice Guideline. The International Guideline. Emily Haesler (Ed). EPUAP/NPIAP/PPPIA: 2019
- (2) Examination of the shear force and the pressure relief capability, Berlincert 2009.
- (3) Pressure readings with Xsensor, 2008.
- (4) CDC (2003), Guidelines for Environmental Infection Control in Health-Care facilities, updated 2019.
- (5) SS-EN ISO 15496:2004, DIN 53.122.
- (6) SS 876 00 20:2017, Healthcare textiles - Mattresses - Specifications and requirements.



Always read the instruction for use prior to usage of the product.

Care of Sweden and CuroCell S.A.M. are trademarks of Care of Sweden AB.
© Care of Sweden AB, 2016. All rights reserved. Last update: 2016-11-15.



Mattress selector

- Guidelines when choosing a product
See how to choose the right mattress at:
www.careofsweden.com/mattress-selector

 care of sweden
Supporting Life

Care of Sweden AB, Box 146, SE-514 23 Tranemo, Sweden
+46 (0)771 106 600 | info@careofsweden.se | www.careofsweden.com