

# CuroCell S.A.M.® PRO

Extraordinary comfort and durability combined



CuroCell S.A.M PRO with cover Steel

# CuroCell S.A.M.® PRO

The CuroCell S.A.M.® PRO overlay mattress can be used as an effective aid in the prevention and in the treatment of pressure ulcers up to and including category 3<sup>(1)</sup>. It can also be used as an aid for pain therapy. For additional comfort, the model CuroCell S.A.M.® PRO CF10 is also available. This model includes an extra soft inner cover that increases the comfort.

This constant low pressure mattress offers high level of comfort as well as effective pressure redistribution <sup>(2)</sup>. With the vastly improved higher maximum weight capacity of 250 kg, this overlay will meet the ever increasing challenges of heavier patients.

The CuroCell S.A.M.® PRO features a new cost effective pump that includes both CPR and Transport functions. The pump automatically determines when

to change pressure in contrast to running constantly. This makes the pump silent running, vibration free, power efficient and long lasting.

Also available is the CuroCell S.A.M.® PRO CF16, which is a mattress replacement system equipped with an integrated foam core, which is placed directly onto the bed frame.

## REFERENCES


- (1) \*Haesler, E. (Ed). National Pressure Ulcer Advisory Panel, European Pressure Ulcer Advisory Panel and Pan Pacific Pressure Injury Alliance. (2014). Prevention and Treatment of Pressure Ulcers: Clinical Guideline. Cambridge Media: Perth, Australia.  
 (2) Examination of the shear force and the pressure relief capability, Berlin cert 2009.  
 (3) ISO 16603 "resistance of penetration by blood and body fluids", ISO 16604 "resistance of penetration by blood-borne pathogens", EN 14126.  
 (4) SS-EN ISO 15496:2004, DIN 53122-1.



Always read the instruction for use prior to usage of the product.

Care of Sweden and CuroCell S.A.M. © are trademarks of Care of Sweden AB. © Care of Sweden AB, 2019. All rights reserved. Last update: 2019-09-11.



FEATURES & CHARACTERISTICS	FACTS
<ul style="list-style-type: none"> <li>• CPR function – rapid deflation</li> <li>• Soft, elastic, vapour permeable<sup>(4)</sup> hygienic cover with welded seams</li> <li>• Very good cleaning properties</li> <li>• Silent running and durable pump</li> <li>• Easy to handle and maintain</li> <li>• Elastic fastening band</li> <li>• Energy efficient</li> <li>• Maximum Pressure (caring mode)</li> <li>• Comfort settings on dial, 0-250 kg</li> </ul>	<ul style="list-style-type: none"> <li>• <b>Type: Constant low pressure mattress</b></li> <li>• <b>Height: 10 cm</b></li> <li>• <b>Function: Prevention and treatment up to and including category 3<sup>(1)</sup></b></li> <li>• <b>Recommended user weight: 0-250 kg</b></li> <li>• <b>Input voltage: 230 V/50 Hz</b></li> <li>• <b>Sound level, pump, max: 28 dB</b></li> <li>• <b>Warranty: 2 years</b></li> <li>• <b>Machine wash max 95 °C, tumble drying</b></li> <li>• <b>Cleaning of cover: Wipe with cleaning agent and/or disinfectants</b></li> </ul>
OPTIONAL	SIZE
<ul style="list-style-type: none"> <li>• Extra soft inner cover (CF10) model</li> <li>• Durable bag for transportation and storage</li> </ul>	<p>80/85/90/105/120x200/210 cm.                      If you need any other sizes please contact your distributor.</p>
 <p><b>SELECT YOUR MATTRESS</b>                      - Visit our website for guidelines when choosing a mattress.  <a href="http://www.careofsweden.com/mattress-selector">www.careofsweden.com/mattress-selector</a></p>	