

CuroCell S.A.M[®] PRO CF16

Versatile & durable mattress replacement system
with high level of comfort



CuroCell S.A.M PRO CF16 with cover Olivia/CC

CuroCell S.A.M.[®] PRO CF16

The CuroCell S.A.M.[®] PRO CF16 mattress replacement system with an integrated foam core can be used as an effective aid in the prevention and in the treatment of pressure ulcers up to and including category 3⁽¹⁾. It can also be used as an aid for pain therapy.

This constant low pressure mattress offers an exceptionally high level of comfort as well as excellent pressure redistribution⁽²⁾. With the vastly improved higher maximum weight capacity of 250 kg, this overlay will meet the ever increasing challenges of heavier patients.

The CuroCell S.A.M.[®] PRO CF16 features a new cost effective pump that includes both CPR and Transport functions. The pump automatically determines when to change pressure in contrast to running constantly.

This does not only make the pump virtually silent but it also makes it durable and long lasting.

Thanks to the extra soft inner cover, the already low pressure⁽³⁾ is further reduced.

REFERENCES

(1) *Haesler, E. (Ed). National Pressure Ulcer Advisory Panel, European Pressure Ulcer Advisory Panel and Pan Pacific Pressure Injury Alliance. (2014). Prevention and Treatment of Pressure Ulcers: Clinical Guideline. Cambridge Media: Perth, Australia.

(2) Examination of the shear force and the pressure relief capability, Berlin Cert 2009.

(3) Pressure readings with Xsensor, 2008.

(4) ISO 16603 "resistance of penetration by blood and body fluids", ISO 16604 "resistance of penetration by blood-borne pathogens", EN 14126.


(5) SS-EN ISO 15496:2004, DIN 53122-1.



Always read the instruction for use prior to usage of the product.

Care of Sweden and CuroCell S.A.M. © are trademarks of Care of Sweden AB. © Care of Sweden AB, 2019. All rights reserved. Last update: 2019-08-08.



FEATURES & CHARACTERISTICS	FACTS
<ul style="list-style-type: none"> • CPR function – rapid deflation • Soft, elastic, vapour permeable⁽⁵⁾ hygienic cover • Very good cleaning properties • Silent running and durable pump • Easy to handle and maintain • Unique labelling • Energy efficient • Maximum Pressure (caring mode) • Comfort settings on dial, 0-250 kg • Extra soft inner cover • Heel function • Durable bag for transportation and storage 	<ul style="list-style-type: none"> • Type: Constant low pressure mattress • Height: 16 cm • Function: Prevention and treatment up to and including category 3⁽¹⁾ • Recommended user weight: 0-250 kg • Input voltage: 230 V/50 Hz • Sound level, pump, max: 28 dB • Warranty: 2 years • Machine wash max 95 °C, tumble drying • Cleaning of cover: Wipe with cleaning agent and/or disinfectants
OPTIONAL	SIZE
<ul style="list-style-type: none"> • Hygiene cover with welded seams 	80/85/90/105/120x200/210 cm. If you need any other sizes please contact your distributor.
 <p>SELECT YOUR MATTRESS - Visit our website for guidelines when choosing a mattress. www.careofsweden.com/mattress-selector</p>	