

CuroCell S.A.M.® Switch



Flexible mattress overlay system
with very high comfort

CuroCell S.A.M.® Switch

The CuroCell S.A.M.® Switch is used as an aid in the prevention and treatment of pressure ulcers up to category 2⁽¹⁾, and for pain therapy. It is a static air mattress overlay system that provides very good comfort and pressure⁽²⁾ redistribution.

CuroCell S.A.M.® Switch has the added advantage that the pump can be used as a support pump only, making the mattress system both noiseless and vibration free.

The mattress is made from very soft and durable materials and comes with an easy to use, very quiet pump. Thanks to the extra soft inner cover (accessory), the already low pressure⁽³⁾ is further reduced.

Also available is the CuroCell S.A.M.® Switch CF16, which is a mattress replacement system equipped with an integrated foam core, and is placed directly onto the bed frame.

Care of Sweden's mattresses are supplied with a four way stretch hygienic outer cover which conforms to relevant standards for hygiene⁽⁴⁾, quality and safety. All covers are manufactured from modern, recyclable materials and are vapour permeable⁽⁵⁾, which reduce the risk of skin maceration.

Facts

Type: Static air mattress
Height: 10 cm
Function: Prevention and treatment up to and including category 2⁽¹⁾
Recommended user weight: 0-200 kg
Power supply: 230 V/50 Hz
Warranty: 2 years
Cleaning of cover: Wipe with cleaning agent and/or disinfectants.
Machine wash max 95 °C, tumble drying.



Features & Characteristics

- Can be used as a support pump only which makes the mattress system both noiseless and vibration free
- Easy to handle and maintain
- Elastic fastening band
- Silent and energy-efficient pump
- Soft, elastic, vapour permeable⁽⁵⁾ hygienic cover
- Unique labelling
- Very good cleaning properties

Sizes

80/85/90/105/120x200/210 cm.
If you need any other sizes please contact your distributor.

Optional

- Bag for the pump
- Durable bag for transportation and storage
- Extra soft inner cover
- Hygiene cover with welded seams
- Transport function

References

- (1) National Pressure Ulcer Advisory Panel, European Pressure Ulcer Advisory Panel and Pan Pacific Pressure Injury Alliance. Prevention and Treatment of Pressure Ulcers: Quick Reference Guide. Emily Haesler (Ed.). Cambridge Media: Osborne Park, Western Australia; 2014.
(2) Examination of the shear force and the pressure relief capability, Berlincert 2009.
(3) Pressure readings with Xsensor, 2008.
(4) ISO 16603 "resistance of penetration by blood and body fluids", ISO 16604 "resistance of penetration by blood-borne pathogens", EN 14126, SS 876 00 19 "Health care textiles – Bacteria penetration – Wet".
(5) SS-EN ISO 15496:2004, DIN 53.122.



Always read the instruction for use prior to usage of the product.

Care of Sweden and CuroCell S.A.M. are trademarks of Care of Sweden AB.
© Care of Sweden AB, 2016. All rights reserved. Last update: 2016-11-15.



Mattress selector

- Guidelines when choosing a product
See how to choose the right mattress at:
www.careofsweden.com/mattress-selector

 care of sweden
Supporting Life

Care of Sweden AB, Box 146, SE-514 23 Tranemo, Sweden
+46 (0)771 106 600 | info@careofsweden.se | www.careofsweden.com