

CuroCell® A4 OP10

A fully automatic and flexible
alternating pressure mattress overlay
system with integrated heel fuction



CuroCell® A4 OP10 with cover Stone

CuroCell® A4 OP10

The CuroCell® A4 OP10 provides simplicity together with user comfort while effectively aiding in the prevention and treatment of pressure ulcers/pressure injuries up to and including category 4, and unclassified pressure ulcers. It is designed for simplicity and provides features for excellent user comfort, reduction of shear forces and improvement of the microclimate.

The CuroCell® A4 OP10 mattress system includes an intelligent pump with a user-friendly interface. It is fully automatic without need for manual settings. The system identifies the user and automatically manages all settings to optimise the pressure according to the user's needs.

All mattresses are supplied with a four way stretch hygienic cover which meets standards for hygiene⁽³⁾, quality and safety. All covers are manufactured from modern, recyclable materials and are vapour permeable⁽⁴⁾, which reduce the risk of skin maceration.



OPTIONAL

- Durable bag
- Sheet holder

SIZES

80/85/90/x/200/210 x 10 cm

REFERENCES

- (1) National Pressure Ulcer Advisory Panel, European Pressure Ulcer Advisory Panel and Pan Pacific Pressure Injury Alliance. Prevention and Treatment of Pressure Ulcers: Quick Reference Guide. Emily Haesler (Ed.). Cambridge Media: Osborne Park, Western Australia; 2014.
- (2) "EN ISO 11201 Acoustics - Noise emitted by machinery and equipment - Determination of emission sound pressure levels", SP 2018.
- (3) ISO 16603 "resistance of penetration by blood and body fluids", ISO 16604 "resistance of penetration by blood-borne pathogens", EN 14126, "performance requirements and tests methods for protective clothing against infective agents".
- (4) SS-EN ISO 15496:2004, DIN 53122-1



Always read the instruction for use prior to usage of the product.

Care of Sweden and CuroCell are trademarks of Care of Sweden AB.
© Care of Sweden AB, 2019. All rights reserved. Last update: 2019-05-28.

FEATURES MATTRESS	FACTS
<ul style="list-style-type: none"> • Heel function • Unique CPR-function with transport function • Unique design for easy handling and usage • Well protected zipper • One-piece mattress 	<ul style="list-style-type: none"> • Type: Alternating pressure mattress overlay system • Height: 10 cm • Recommended user weight: 0-160 kg • Input voltage: 100-240 V / 50-60 Hz • Sound level, pump, max.: 17 dBA⁽²⁾ • Warranty: 2 years • PVC free • Cleaning of cover: Wipe with cleaning agent and/or disinfectants • Machine wash max 95 °C, tumble drying • Chlorine resistant (1%) multi stretch cover (cover Stone with welded seams, 10%) <div data-bbox="879 1816 986 1921" style="text-align: center;"> </div> <p data-bbox="1007 1850 1235 1906">Scan QR code to watch presentation video.</p> <div data-bbox="879 1957 986 2063" style="text-align: center;"> </div> <p data-bbox="1007 1951 1362 1980">SELECT YOUR MATTRESS</p> <p data-bbox="1007 1989 1362 2040">- Visit our website for guidelines when choosing a mattress.</p> <p data-bbox="1007 2047 1410 2069">www.careofsweden.com/mattress-selector</p>
<p>FEATURES PUMP</p>	
<ul style="list-style-type: none"> • Fully automated pump that adapts individually • User friendly interface • Audible and visual notifications • Three program modes: Pulsating, Alternating and CLP (constant low pressure) • Caring mode (max-firm) • Pack&Go® function • Silent running and vibration free • Low voltage (12V power supply) • Soft friction hangers with extra soft tip • Air flow control; Constant temperature, increased user comfort, energy reduction 	