

CuroCell® Nova


Simple & functional mattress overlay system



CuroCell Nova with cover Steel

CuroCell® Nova

The CuroCell® Nova is used as an aid in the prevention and treatment of pressure ulcers up to category 2⁽¹⁾. It is a mattress overlay system that is ideal for community and home care and is easy to handle and operate. The one-piece mattress has air cells of polyurethane and a built-in heel function.



REFERENSER


- (1) *Haesler, E. (Ed). National Pressure Ulcer Advisory Panel, European Pressure Ulcer Advisory Panel and Pan Pacific Pressure Injury Alliance. (2014). Prevention and Treatment of Pressure Ulcers: Clinical Guideline. Cambridge Media: Perth, Australia.
- (2) ISO 16603 "resistance of penetration by blood and body fluids", ISO 16604 "resistance of penetration by blood-borne pathogens", EN 14126.
- (3) SS-EN ISO 15496:2004, DIN 53122-1.



Always read the instruction for use prior to usage of the product.

Care of Sweden and CuroCell S.A.M. © are trademarks of Care of Sweden AB.
© Care of Sweden AB, 2019. All rights reserved. Last update: 2019-09-11.



FEATURES & CHARACTERISTICS	FACTS
<ul style="list-style-type: none"> • CPR function – rapid deflation • Soft, elastic, vapour permeable⁽⁴⁾ hygienic cover • Very good cleaning properties • Silent running and durable pump • Easy to handle and maintain • Elastic fastening band • Unique labelling • Energy efficient • Static function and Max firm function • Comfort settings on dial, 0-150 kg 	<ul style="list-style-type: none"> • Type: Alternating pressure mattress replacement system • Height: 13 cm • Function: Prevention and treatment up to and including category 2⁽¹⁾ • Recommended user weight: 0-150 kg • Input voltage: 230 V/50 Hz • Sound level, pump, max: 21 dBA • Warranty: 2 years • Machine wash max 95 °C, tumble drying • Cleaning of cover: Wipe with cleaning agent and/or disinfectants
OPTIONAL	SIZE
<ul style="list-style-type: none"> • Extra soft inner cover • Hygienic cover with welded seams • Durable bag for transportation and storage 	<p>80x200 cm.</p>
	 <p>SELECT YOUR MATTRESS - Visit our website for guidelines when choosing a mattress. www.careofsweden.com/mattress-selector</p>