

CuroCell® Nova CX13



Dynamic pressure care mattress overlay system

CuroCell® Nova CX13

The CuroCell® Nova CX13 is used as an aid in the prevention and treatment of pressure ulcers up to category 2⁽¹⁾ and is ideal for community and home care. It is a mattress overlay system that is easy to handle and operate.

The CuroCell® Nova CX13 uses the alternating pressure method and is available in several different models. It has replaceable polyurethane cells, which are housed in a separate inner cover.

Also available is CuroCell® Nova CX19, a mattress replacement system equipped with an integrated foam core.

Care of Sweden's mattresses are supplied with a four way stretch hygienic outer cover which conforms to relevant standards for hygiene⁽²⁾, quality and safety. All covers are manufactured from modern, recyclable materials and are vapour permeable⁽³⁾, which reduce the risk of skin maceration.

Facts

Type: Alternating pressure mattress

Height: 13 cm

Function: Prevention and treatment up to and including category 2⁽¹⁾

Recommended user weight: 0-150 kg

Power supply: 230 V/50 Hz

Noise level of pump: 21 dBA

Warranty: 2 years

Cleaning of cover: Wipe with cleaning agent and/or disinfectants.

Machine wash max 95 °C, tumble drying.



Features & Characteristics

- CPR - Rapid deflation
- Easy handling and cleaning
- Replaceable cells in polyurethane
- Silent, energy-saving pump
- Soft, elastic, vapour permeable⁽³⁾ hygienic cover
- Transport function

Sizes

80/85/90/105/120x200 cm and 80/85/90x210 cm.

If you need any other sizes please contact your distributor.

Optional

- Durable bag
- Heel function
- Hygienic cover with welded seams
- Ventilated cells (Low Air Loss)

References

(1) National Pressure Ulcer Advisory Panel, European Pressure Ulcer Advisory Panel and Pan Pacific Pressure Injury Alliance. Prevention and Treatment of Pressure Ulcers: Quick Reference Guide. Emily Haesler (Ed.). Cambridge Media: Osborne Park, Western Australia; 2014.

(2) ISO 16603 "resistance of penetration by blood and body fluids", ISO 16604 "resistance of penetration by blood-borne pathogens", EN 14126, SS 876 00 19 "Health care textiles – Bacteria penetration – Wet".

(3) SS-EN ISO 15496:2004, DIN 53.122.



Always read the instruction for use prior to usage of the product.

Care of Sweden and CuroCell are trademarks of Care of Sweden AB.

© Care of Sweden AB, 2016. All rights reserved. Last update: 2016-11-22.



Mattress selector

- Guidelines when choosing a product
See how to choose the right mattress at:
www.careofsweden.com/mattress-selector

 care of sweden
Supporting Life

Care of Sweden AB, Box 146, SE-514 23 Tranemo, Sweden
+46 (0)771 106 600 | info@careofsweden.se | www.careofsweden.com