

CuroCell® Nova CX13

Dynamic pressure care mattress overlay system



CuroCell® Nova CX13 with cover Olivia

CuroCell® Nova CX13

The CuroCell® Nova CX13 is used as an aid in the prevention and treatment of pressure ulcers up to category 2⁽¹⁾ and is ideal for community and home care. It is a mattress overlay system that is easy to handle and operate.

The CuroCell® Nova CX13 uses the alternating pressure method and is available in several different models. It has replaceable polyurethane cells, which are housed in a separate inner cover.

Also available is CuroCell® Nova CX19, a mattress replacement system equipped with an integrated foam core.

Care of Sweden's mattresses are supplied with a four way stretch hygienic outer cover which conforms to relevant standards for hygiene⁽²⁾, quality and safety.

All covers are manufactured from modern, recyclable materials and are vapour permeable⁽³⁾, which reduce the risk of skin maceration.

REFERENSER

(1) *Haesler, E. (Ed). National Pressure Ulcer Advisory Panel, European Pressure Ulcer Advisory Panel and Pan Pacific Pressure Injury Alliance. (2014). Prevention and Treatment of Pressure Ulcers: Clinical Guideline. Cambridge Media: Perth, Australia.

(2) ISO 16603 "resistance of penetration by blood and body fluids", ISO 16604 "resistance of penetration by blood-borne pathogens", EN 14126.


(3) SS-EN ISO 15496:2004, DIN 53122-1.



Always read the instruction for use prior to usage of the product.

Care of Sweden and CuroCell are trademarks of Care of Sweden AB.
© Care of Sweden AB, 2020. All rights reserved. Last update: 2020-02-04.



FEATURES & CHARACTERISTICS	FACTS
<ul style="list-style-type: none"> • CPR function – rapid deflation • Easy handling and cleaning • Replaceable cells in polyurethane • Silent, energy-saving pump • Soft, elastic, vapour permeable⁽³⁾ hygienic cover • Transport function 	<ul style="list-style-type: none"> • Type: Alternating pressure mattress overlay • Height: 13 cm • Function: Prevention and treatment up to and including category 2⁽¹⁾ • Recommended user weight: 0-150 kg • Input voltage: 230 V/50 Hz • Sound level, pump, max: 21 dBA • Warranty: 2 years • Machine wash max 95 °C, tumble drying • Cleaning of cover: Wipe with cleaning agent and/or disinfectants
<p>OPTIONAL</p> <ul style="list-style-type: none"> • Durable bag • Heel function • Hygienic cover with welded seams 	
 <p>SELECT YOUR MATTRESS - Visit our website for guidelines when choosing a mattress. www.careofsweden.com/mattress-selector</p>	<p>SIZES</p> <p>80/85/90/105/120x200 cm and 80/85/90x210 cm. If you need any other sizes please contact your distributor.</p>