

CuroCell® Nova CX19

Dynamic pressure care mattress replacement system



CuroCell® Nova CX19 with cover Olivia

CuroCell® Nova CX19

CuroCell® Nova CX19 is used as an aid in the prevention and treatment of pressure ulcers up to category 2⁽¹⁾. It is used in all types of healthcare environments and is a mattress replacement system that is very easy to handle and operate.

CuroCell® Nova CX19 works according to the alternating pressure principle and are available in several different versions. It has replaceable cells in polyurethane which are placed in a separate inner cover and is equipped with an integrated foam core.

Care of Sweden's mattresses are supplied with a four way stretch hygienic outer cover which conforms to relevant standards for hygiene⁽²⁾, quality and safety. All covers are manufactured from modern, recyclable materials and are vapour permeable⁽³⁾, which reduce the risk of skin maceration.

REFERENSER


- (1) *Haesler, E. (Ed). National Pressure Ulcer Advisory Panel, European Pressure Ulcer Advisory Panel and Pan Pacific Pressure Injury Alliance. (2014). Prevention and Treatment of Pressure Ulcers: Clinical Guideline. Cambridge Media: Perth, Australia.
 (2) ISO 16603 "resistance of penetration by blood and body fluids", ISO 16604 "resistance of penetration by blood-borne pathogens", EN 14126.
 (3) SS-EN ISO 15496:2004, DIN 53122-1.



Always read the instruction for use prior to usage of the product.

Care of Sweden and CuroCell are trademarks of Care of Sweden AB.
 © Care of Sweden AB, 2020. All rights reserved. Last update: 2020-02-04.



FEATURES & CHARACTERISTICS	FACTS
<ul style="list-style-type: none"> • CPR function – rapid deflation • Easy handling and cleaning • Heel function • Replaceable cells in polyurethane • Silent, energy-saving pump • Soft, elastic, vapour permeable⁽³⁾ hygienic cover • Transport function 	<ul style="list-style-type: none"> • Type: Alternating pressure mattress replacement system • Height: 19 cm • Function: Prevention and treatment up to and including category 2⁽¹⁾ • Recommended user weight: 0-150 kg • Fire: EN 597-1, -2 • Input voltage: 230 V/50 Hz • Sound level, pump, max: 21 dBA • Warranty: 2 years • Machine wash max 95 °C, tumble drying • Cleaning of cover: Wipe with cleaning agent and/or disinfectants
<p>OPTIONAL</p> <ul style="list-style-type: none"> • Cable holder • Comfort fabric • Durable bag • Hygienic cover with welded seams • Sheet holder 	
 <p>SELECT YOUR MATTRESS - Visit our website for guidelines when choosing a mattress. www.careofsweden.com/mattress-selector</p>	<p>SIZES</p> <p>80/85/90/105/120x200 cm and 80/85/90x210 cm. If you need any other sizes please contact your distributor.</p>