

# CuroCell® 3 CX13



Effective & flexible mattress overlay system

# CuroCell<sup>®</sup> 3 CX13

CuroCell<sup>®</sup> 3 CX13 is used as an aid in the prevention and treatment of pressure ulcers up to category 3<sup>(1)</sup>. It is used in all types of healthcare environments and is an effective and functional mattress overlay system.

CuroCell<sup>®</sup> 3 CX13 works according to the alternating pressure principle and has replaceable cells in polyurethane which are placed in a separate inner cover.

Also available is CuroCell 3 CX19, a mattress replacement system equipped with an integrated foam core.

Care of Sweden's mattresses are supplied with a four way stretch hygienic outer cover which conforms to relevant standards for hygiene, quality and safety<sup>(2)(4)</sup>. All covers are manufactured from modern, recyclable materials and are vapour permeable<sup>(3)</sup>, which reduce the risk of skin maceration.

## Facts

Type: Alternating pressure mattress

Height: 13 cm

Function: as an aid for prevention and treatment up to and including category 3<sup>(1)</sup>

Recommended user weight: up to 160 kg

Power supply: 230 V/50 Hz

Noise level of pump: 21 dBA

Warranty: 2 years

Cleaning of cover: Wipe with cleaning agent and/or disinfectants.

Machine wash max 95 °C, tumble drying.



## Features & Characteristics

- CPR function – rapid deflation
- Digital pump with built-in sensor
- Replaceable air cells of polyurethane
- Silent and energy-saving pump
- Soft, elastic, vapour permeable<sup>(3)</sup> hygienic cover
- Sound and light notification for low pressure
- Static function with automatic return
- Transport function
- Good cleaning properties

## Sizes

80/85/90/105/120x200 cm and 80/85/90x210 cm.

If you need any other sizes please contact your distributor.

## Optional

- Durable bag
- Heel function
- Hygienic cover with welded seams

## References

(1) European Pressure Ulcer Advisory Panel, National Pressure Injury Advisory Panel and Pan Pacific Pressure Injury Alliance. Prevention and Treatment of Pressure Ulcers/Injuries: Clinical Practice Guideline. The International Guideline. Emily Haesler (Ed). EPUAP/NPIAP/PPPIA: 2019

(2) CDC (2003), Guidelines for Environmental Infection Control in Health-Care facilities, updated 2019.

(3) SS-EN ISO 15496:2004, DIN 53.122.

(4) SS 876 00 20:2017, Healthcare textiles - Mattresses - Specifications and requirements.



Always read the instruction for use prior to usage of the product.

Care of Sweden and CuroCell are trademarks of Care of Sweden AB.

© Care of Sweden AB, 2021. All rights reserved. Last update: 2021-06-03.



### Mattress selector

- Guidelines when choosing a product  
See how to choose the right mattress at:  
[www.careofsweden.com/mattress-selector](http://www.careofsweden.com/mattress-selector)

 care of sweden  
Supporting Life

Care of Sweden AB, Box 146, SE-514 23 Tranemo, Sweden  
+46 (0)771 106 600 | [info@careofsweden.se](mailto:info@careofsweden.se) | [www.careofsweden.com](http://www.careofsweden.com)