


CuroCell® 4 AD

Comfortable alternating pressure mattress system



CuroCell 4 AD with Top part Olivia/Bottom part CC

CuroCell® 4 AD


The CuroCell® 4 AD is used as an aid in the prevention and treatment of pressure ulcers up to and including category 4⁽¹⁾. It's a highly functional, comfortable, cost-effective mattress replacement system and is ideal for all healthcare environments.

Care of Sweden's mattresses are supplied with a four way stretch hygienic outer cover which conforms to relevant standards for hygiene⁽²⁾, quality and safety.

All covers are manufactured from modern, recyclable materials and are vapour permeable⁽³⁾, which reduce the risk of skin maceration.

THE ACCESSORY PACKAGE PRO ALSO INCLUDES THE OPTIONS

- Heel function
- Cable holder
- Comfort Fabric
- Sheet holder
- Durable bag



(Accessory)

REFERENCES


- (1) *Haesler, E. (Ed). National Pressure Ulcer Advisory Panel, European Pressure Ulcer Advisory Panel and Pan Pacific Pressure Injury Alliance. (2014). Prevention and Treatment of Pressure Ulcers: Clinical Guideline. Cambridge Media: Perth, Australia.
- (2) ISO 16603 "resistance of penetration by blood and body fluids", ISO 16604 "resistance of penetration by blood-borne pathogens", EN 14126
- (3) SS-EN ISO 15496:2004, DIN 53122-1.



Always read the instruction for use prior to usage of the product.



Care of Sweden and CuroCell are trademarks of Care of Sweden AB.
© Care of Sweden AB, 2019. All rights reserved.
Last update: 2019-09-11.

FEATURES & CHARACTERISTICS	FACTS
<ul style="list-style-type: none"> • CPR function – rapid deflation • Digital pump with built-in sensor • Easy handling • Handle for transport and moving • Replaceable air cells of polyurethane with Cell-in-Cell system • Silent and energy-saving pump • Soft, elastic, vapour permeable⁽³⁾ hygienic cover • Sound and light alarm for low pressure • Sound and light alarm for power failure • Static function with automatic return • Transport function • Very good cleaning properties 	<ul style="list-style-type: none"> • Type: Alternating pressure mattress • Height: 21 cm • Function: Prevention and treatment up to and including category 4⁽¹⁾ • Recommended user weight: 0-200 kg • Power supply: 230 V/50 Hz • Noise level of pump: 27 dBA • Warranty: 2 years • Cleaning of cover: Wipe with cleaning agent and/or disinfectants. • Machine wash max 95 °C, tumble drying.
OPTIONAL	SIZES
<ul style="list-style-type: none"> • Cable holder (CuroCell 4® PRO) • Comfort fabric (CuroCell 4® PRO) • Durable bag (CuroCell 4® PRO) • Heel function (CuroCell 4® PRO) • Sheet holder (CuroCell 4® PRO) • Hygienic cover with welded seams • Low friction cover • Ventilated cells (Low Air Loss) 	<p>80/85/90/105/120x200 cm and 80/85/90x210 cm. If you need any other sizes please contact your distributor.</p>
	 <p>MATTRESS SELECTOR – Guidelines when choosing a product See how to choose the right mattress at: www.careofsweden.com/mattress-selector</p>