

CuroCell Cirrus® 2.0

Pain therapy pressure care system



CuroCell Cirrus® 2.0 with Top part Olivia/ bottom part CC

CuroCell Cirrus® 2.0

CuroCell Cirrus® 2.0 is used as an aid in the prevention and treatment of pressure ulcers up to category 4 (0-230 kg)⁽¹⁾, and for pain therapy. It is a replacement system ideal for all healthcare environments.

CuroCell Cirrus® 2.0 combines air with a soft foam surface and provides high user comfort and pressure redistribution⁽²⁾. The unique wave like alternating process is particularly suitable for pain therapy. The added heel function offers an even better protection by relieving pressure on the heels.

The silent and easy to use CuroCell Cirrus® pump with four different settings; alternating, constant low pressure, alternating + seat inflation and constant low pressure + seat inflation make this mattress very functional for all patient groups. In addition the dual cover system reduces the shearing effect when adjusting the bed.

Care of Sweden's mattresses are supplied with a four way stretch hygienic outer cover which con-

forms to relevant standards for hygiene⁽³⁾, quality and safety. All covers are manufactured from modern, recyclable materials and are vapour permeable⁽⁴⁾, which reduce the risk of skin maceration.

OPTIONAL

- Durable bag
- Low friction cover

REFERENCES

(1) *Haesler, E. (Ed). National Pressure Ulcer Advisory Panel, European Pressure Ulcer Advisory Panel and Pan Pacific Pressure Injury Alliance. (2014). Prevention and Treatment of Pressure Ulcers: Clinical Guideline. Cambridge Media: Perth, Australia.

(2) Examination of the shear force and the pressure relief capability, Berlin Cert 2009.

(3) ISO 16603 "resistance of penetration by blood and body fluids", ISO 16604 "resistance of penetration by blood-borne pathogens", EN 14126.

(4) SS-EN ISO 15496.



Always read the instruction for use prior to usage of the product.

Care of Sweden, CuroCell Cirrus are trademarks of Care of Sweden AB. © Care of Sweden AB, 2019. All rights reserved.
Last update: 2019-08-29.



FEATURES & CHARACTERISTICS	FACTS
<ul style="list-style-type: none"> • Alternating pressure on soft foam surface • Cabel holder • Digital and easy to use pump • Easy handling and cleaning • Heel slope, for pressure relief of heels • Hygienic cover with welded seams • Rapid deflation of the mattress • Silent and energy-saving pump • Soft, elastic, vapour permeable⁽⁴⁾ hygienic cover • Sound and light notification for low pressure • Sound and light notification for power failure • CLP-function with automatic return • Transport function 	<ul style="list-style-type: none"> • Type: Alternating pressure mattress • Height: 21 cm • Function: Prevention and treatment up to and including category 4⁽¹⁾ • Recommended user weight: 0-230 kg • Fire: EN 597-1, -2 • Power supply: 230 V/50 Hz • Noise level of pump: 27 dBA • Warranty: 2 years • Cleaning of cover: Wipe with cleaning agent and/or disinfectants • Machine wash max 95 °C, tumble drying <div data-bbox="884 1816 991 1921"> </div> <p data-bbox="1011 1845 1238 1899">Scan QR code to watch presentation video.</p>
<p>SIZE</p>	<div data-bbox="884 1951 991 2056"> </div> <p>SELECT YOUR MATTRESS</p> <p>- Visit our website for guidelines when choosing a mattress. www.careofsweden.com/mattress-selector</p>
<p>80/85/90/105/120x200 cm and 80/85/90x210 cm. If you need any other sizes please contact your distributor.</p>	