

Curera[®] Neck Pillow

Comfortable neck relief



 **curera**[™]
positionering

 care of sweden

SUPPORTING LIFE

Curera® Neck Pillow

If a person makes too few spontaneous movements, there is an increased risk of pressure ulcers. Regular changes of position reduce the risk of pressure ulcers. Position as required in both reclining and seated position.

Curera® Positioning is a range of stable pillows, rolls, and wedges that meet most needs for support and relief. They are also a great help when changing dressings and for patient care.

Our soft hygiene covers and cotton covers are also available as accessories to the pillows and rolls.

The hygiene cover does not allow liquid or air to penetrate, but is vapour permeable (vapour is wicked away).


Make sure that the patient is positioned correctly so that airways are unobstructed.

Working with Curera® Positioning is simple, and the benefits are many. Contact Care of Sweden for training and further information.

POSITIONING EXAMPLE



PRESSURE RELIEF, SUPPORT AND POSITIONING OF HEAD/NECK.

CLEANING	FACTS
<p>The Hygiene Cover can be wiped clean with detergent and/or disinfectant. Machine wash at max. 95 °C.</p> <p>Inner pillow: Machine wash at max. 95 °C. Can be tumble-dried or dried in an airing cupboard; avoid over-drying. For faster drying, run the spin cycle three or four times.</p> <p>Cotton cover: Machine wash at max. 70 °C.</p>	<ul style="list-style-type: none"> • Filling: 50 % HR foam/ 50 % polyester fibre • Warranty: 1 year • Sizes: 32, 38, 50 cm • Accessories: Hygiene Cover: Polyurethane/polyester Cotton cover <p> Read the user manual prior to using the product.</p>

Care of Sweden and Curera are registered trademarks of Care of Sweden AB. © Care of Sweden AB, 2019. All rights reserved. Reproduction prohibited. Last updated: 2019-07-03