


Curera® Support Wedge



 **curera**™
positioning

For positioning & position changes

Curera® Support Wedge


Curera® Support Wedge can be used to place a patient/user in side position in order to relieve sacrum, which is an extra sensitive area. For the side position, two support wedges are used together.

The support wedge is easy to use and clean. It is delivered with a liquid resistant cover that can be wiped off with an cleaning agent and/or disinfectants. It can also be machine washed.

Support wedges are placed differently depending on the mattress that is used, see picture below or user manual for more information.



Example of positionings



Facts

Filling: 100 % foam
Cover: Polyurethane/polyester
Warranty: 1 year
Cleaning of cover: Wipe with cleaning agent and/or disinfectants.
Machine wash max 95 °C, tumble drying.

Size

70x26x15/5 cm

Example of application

- As kick protection against bed rails
- For relief in prevention of pressure ulcers
- When the user has problems to lie completely on the side
- When the user has problems with position changes
- For other positioning



Always read the instruction for use prior to usage of the product.

Care of Sweden and Curera are trademarks of Care of Sweden AB.
© Care of Sweden AB, 2017. All rights reserved. Last update: 2017-02-28.



Scan QR code to watch instructional video.