

Optimal Evacuation Sheet



For quick & simple evacuation

Optimal Evacuation Sheet

Optimal Evacuation Sheet makes it possible to quickly and simply move bedridden persons away from danger. It can also be used for evacuation in stairs.

It is placed under the existing mattress and is fastened to the mattress with elastic straps. It is easy to wrap the person up in the bed and is available in several sizes in order to fit most mattresses.

Facts

Material: Polyamide/polyester

Fire: EN 597-1, -2

Warranty: 2 years

Cleaning: Wipe with cleaning agent and/or disinfectants. Machine wash max 95 °C, tumble drying.



Sizes

Mattress width: 80, 90, 105 & 120 cm

Mattress length: 200 cm

Optimal Evacuation Sheet is ergonomically designed and produced from very strong material that slides well on both wet and dry surfaces. It is X-ray transparent and equipped with handles for sideways movement. User instructions are included.

A simple and safe means for evacuation. Fulfils requirements of the SS 8760039 (Swedish Standard) for rescue sheets.

Features & Characteristics

- Corner elastic strap
- Easy to clean and handle
- Extensible velcro fastening
- Pull handles on short sides
- Side handles



Always read the instruction for use prior to usage of the product.

Care of Sweden is a trademark of Care of Sweden AB.

© Care of Sweden AB, 2014. All rights reserved. Last update: 2014-08-01.



Scan QR code to watch instructional video.