

# Optimal Hygiene Pillow

Comfortable & durable pillows



# Optimal Hygiene Pillow

For good comfort, we can offer a wide range of accessories, including pillows and pillow covers.

The filling in Optimal Hygiene Pillow consists of ball fibre that can be easily fluffed up to maintain the original shape of the pillow. Thanks to the included hygiene cover, which is removable, the pillow is easy to use and maintain.

Optimal Sens is an anatomical sleeping pillow. It is made of a special foam that conforms to the shape of the neck and head, which provides good support and relief to the neck and head .

A removable hygiene cover is included.

Optimal Pillow Case is durable and long-lasting. Hygiene cover which conforms to relevant hygienic standards<sup>(1)(3)</sup>. Removable with zip.



The hygiene cover does not allow liquid or air to penetrate, but is vapour permeable<sup>(2)</sup> (vapour is wick-ed away). Make sure that the patient is positioned correctly so that airways are unobstructed. The cover can be easily cleaned by wiping or machine washing, which minimises the need to wash the entire pillow.

## REFERENCES

(1) CDC (2003), Guidelines for Environmental Infection Control in Health-Care facilities, updated 2019

(2) SS-EN ISO 15496:2004, DIN 53.122.

(3) SS 876 00 20:2017, Healthcare textiles - Mattresses - Specifications and requirements

OPTIONAL EXTRAS	FACTS
<ul style="list-style-type: none"> <li>• Optimal Hygienic Pillow &amp; Pillow Case</li> <li>• Plastic zip</li> <li>• Welded seams</li> </ul>	<ul style="list-style-type: none"> <li>• <b>Optimal Hygienic Pillow &amp; Pillow Case</b></li> <li>• <b>Material: Polyester/cotton</b></li> <li>• <b>Pillow case: Polyurethane/polyester</b></li> <li>• <b>Cleaning: Machine wash at max. 60 °C (pillow)</b></li> <li>• <b>Optimal Sens Ergonomic Pillow &amp; Pillow Case</b></li> <li>• <b>Core: Viscose elastic foam, 50 kg/m<sup>3</sup></b></li> <li>• <b>Pillow case: Polyurethane/polyester</b></li> <li>• <b>Warranty: 1 year</b></li> <li>• <b>Cleaning the cover: Wipe clean with detergent and/or disinfectant.</b></li> <li>• <b>Machine wash at max. 95 °C.</b></li> </ul> <div style="text-align: center; margin-top: 20px;">  </div> <div style="margin-top: 10px;">  <span>Read the instructions for use prior to using the product.</span> </div>
<h3>SIZES</h3>	
<p><b>Optimal Hygienic Pillow &amp; Pillow Case</b> 35x45 cm 50x60 cm 50x70 cm</p> <p><b>Optimal Sens &amp; Pillow Case</b> 30x48 cm</p> <p><b>Optimal Pillow Case Cover</b> 35x45 cm 50x60 cm 50x70 cm</p>	

Care of Sweden and Curera are registered trademarks of Care of Sweden AB. © Care of Sweden AB, 2021. All rights reserved. Reproduction prohibited. Last updated: 2021-05-31