

# Optimal Mattress Cover



A durable mattress cover  
with a good fit & flexibility

# Optimal Mattress Cover

For added comfort in every unique situation, we can offer a wide range of accessories. Our Optimal Mattress Cover is developed for the special needs in healthcare environments, but is also suitable in a home environment.

Optimal Mattress Cover is an easy, strong and functional hygiene cover with a good fit and flexibility.

## Facts

---

Warranty: 2 years

Cleaning: Wipe with cleaning agent and/or disinfectants. Machine wash max 95 °C, tumble drying.



## Features & Characteristics

---

- Easy to clean
- Elastic on the bottom
- Hygienic
- Liquid resistant
- Many different sizes
- Soft and elastic
- Strong
- Vapour permeable

Optimal Mattress Cover is detachable and has a very well documented hygienic standard<sup>(1)</sup>. The soft and elastic material is vapour permeable<sup>(2)</sup>, which reduces the risk of skin maceration.

## Sizes

---

80x205 cm

85x200 cm

90x200 cm

105x200 cm

120x200 cm

## References

---

(1) ISO 16603 "resistance of penetration by blood and body fluids", ISO 16604 "resistance of penetration by blood-borne pathogens", EN 14126, SS 876 00 19 "Health care textiles – Bacteria penetration – Wet".

(2) SS-EN ISO 15496:2004, DIN 53.122.



Always read the instruction for use prior to usage of the product.

Care of Sweden is a trademark of Care of Sweden AB.

© Care of Sweden AB, 2013. All rights reserved. Last update: 2013-11-27.