

# Optimal Mattress Cover

Durable mattress topper with good fit & flexibility



# Optimal Mattress Cover

For good comfort we offer a wide range of accessories, including the Optimal Mattress Cover, which is designed to meet healthcare needs but can also be used at home.

Optimal Mattress Cover is a simple, durable and functional mattress topper with a good fit and flexibility. Optimal Mattress Cover is removable and has a good hygienic standard<sup>(1)(3)</sup>. The elastic cover does not allow liquid or air to penetrate and is vapour permeable<sup>(2)</sup> (vapour is wicked away), which reduces the risk of skin maceration.

## REFERENCES

(1) CDC (2003), Guidelines for Environmental Infection Control in Health-Care facilities, updated 2019

(2) SS-EN ISO 15496:2004, DIN 53.122.

(3) SS 876 00 20:2017, Healthcare textiles - Mattresses - Specifications and requirements

FUNCTIONS AND PROPERTIES	FACTS
<ul style="list-style-type: none"><li>• Good fit</li><li>• Hygienic</li><li>• Easy to clean</li><li>• Elastic on top and underside</li><li>• Assorted sizes</li><li>• Durable</li><li>• Liquid-resistant</li><li>• Vapour permeable</li></ul>	<ul style="list-style-type: none"><li>• <b>Warranty: 1 year</b></li><li>• <b>Hygiene topper: Polyurethane/polyester</b></li><li>• <b>Cleaning: Wipe clean with detergent and/or disinfectant. Machine wash max. 95°C, tumble dry.</b></li></ul> <div data-bbox="874 1339 960 1406"></div> <div data-bbox="874 1440 933 1496"></div> <p data-bbox="946 1451 1404 1478">Read the instructions for use prior to using the product.</p>
SIZES	
80x205 cm 85x200 cm 90x200 cm 105x200 cm 120x200 cm	

Care of Sweden and Curera are registered trademarks of Care of Sweden AB. © Care of Sweden AB, 2021. All rights reserved. Reproduction prohibited. Last updated: 2021-06-09