

# Curera<sup>®</sup> Ear Pillow

For support and relief of sensitive ears



Ear pillow with cotton cover/hygiene cover

# Curera<sup>®</sup> Ear Pillow

A pillow designed to support and relieve sensitive ears, an area that is particularly prone to pressure ulcers when a user lies still in bed.

The pillow supports the face and the oval-shaped hole totally relieves pressure off the ear.

The ear pillow is sewn in a hygienic fabric that does not allow liquid or air to penetrate, but is vapour permeable (vapour is transported away). Make sure that the patient is positioned correctly so that airways are unobstructed.

The ear pillow comes with a cotton cover. Accordingly, it ought to be used as a product intended for a single user.



Always read the instruction for use prior to usage of the product.

Care of Sweden and Curera are registered trademarks of Care of Sweden AB. © Care of Sweden AB, 2018. All rights reserved. Reproduction prohibited. Last updated: 2023-02-27.



CLEANING	FACTS
<p><b>Pillow</b> The pillow can be wiped clean with detergent and/or disinfectant. Machine wash at max. 95°C.</p> <p><b>Cotton cover</b> Machine wash at max. 70°C.</p>	<ul style="list-style-type: none"><li>• <b>Filling: 50% HR foam/50% polyester fibre</b></li><li>• <b>Warranty: 1 year</b></li><li>• <b>Size: 32x32 cm</b></li></ul>