



CuroCell S.A.M.® PRO

Good comfort, durability, and easy to use

 care of sweden
SUPPORTING LIFE

Good comfort and ease of use with CuroCell S.A.M.® PRO

Thanks to our continuous development efforts, conducted in collaboration with healthcare professionals and users, we have developed a system for the CuroCell S.A.M.® -series based on the need to equalise pressure in different situations: low noise level, user comfort and ease of use.



CuroCell S.A.M.® PRO CF10

Good user comfort

CuroCell S.A.M.® PRO offers a calm environment for users, with the comfort dial enabling pressure distribution adapted to the individual. The system offers low noise- and vibration levels that contributes to a good user comfort. The system is used as an aid to prevent and treat pressure ulcers up to category III.

- **SOFT INNER COVER (OPTIONAL)**
- **GOOD COMFORT**
- **LOW NOISE LEVELS**

CuroCell S.A.M.® PRO is available in three versions

CuroCell S.A.M.® PRO	10 cm overlay mattress up to 250 kg
CuroCell S.A.M.® PRO CF10	10 cm overlay mattress with inner cover up to 250 kg
CuroCell S.A.M.® PRO CF16	16 cm replacement mattress with inner cover up to 250 kg

Digital pump for good durability

Mattress systems used on a daily basis for extended periods and in different situations must be both functional and reliable. This is why CuroCell SAM® PRO is developed with a pump that is activated only when needed, such as when the user change position. This ensures good durability and reliable operation.

- DURABLE
- USER-FRIENDLY INTERFACE
- SIMPLY SET THE USER'S WEIGHT
- ACTIVE ONLY WHEN NEEDED




CuroCell S.A.M.® PRO CF16



Easy to use

Daily activities involving different types of medical devices place strict requirements on a system's ease of use and hygiene properties. CuroCell S.A.M.® PRO is designed to be easily installed, used, and cleaned. The system features a nursing mode that aims to simplify care and makes moving the user in and out of bed easier. When transporting the user, switching to transport mode detaches the pump from the mattress, which maintains its constant low pressure.

- AUDIO-VISUAL NOTIFICATIONS
- DESIGNED FOR EASE OF USE
- SIMPLIFIES CARE WITH THE NURSING MODE
- EASY TO CLEAN AND TO MAINTAIN GOOD HYGIENE



*“Soft and comfortable
with a silky feel”*

CuroCell S.A.M.® exceeded our expectations.
The mattress is soft and comfortable with a silky feel.
It's also easy for personnel to handle, in terms of both
weight and use. One user who got to use the mattress
exclaimed: It's heavenly, so comfortable!

Nurse responsible for bed management

“Distributes pressure evenly”

I've had high-level tetraplegia since 1975 and have found
many other mattresses hard. With the CuroCell S.A.M.®
PRO, I could finally lower the pressure so that I sank
into the mattress with more evenly distributed pressure.
I'd recommend it to anyone with similar problems!

User

“Alleviated the users pain”

We have a user with severe epilepsy, severe
learning difficulties and increasingly limited mobility.
We upgraded their mattress to a CuroCell S.A.M.®
and this alleviated their pain. After a brief period of
redness, the user's pressure ulcers also cleared up.

Nurse, service home, Arvika Municipality

 care of sweden
SUPPORTING LIFE

Care of Sweden AB, Box 146, SE-514 23 Tranemo, Sweden
+46 771 106 600 | export@careofsweden.se | www.careofsweden.com
ID: SAMP-CE-090-PB-EN