

CuroCell Cirrus® 2.0

A mattress replacement system
with foam and air combined



CuroCell Cirrus® 2.0 with Top part Stone/ bottom part CC

 care of sweden
SUPPORTING LIFE

CuroCell Cirrus® 2.0

CuroCell Cirrus® 2.0 is used as an aid in the prevention and treatment of pressure ulcers up to category 4⁽¹⁾. It is a replacement system intended for all healthcare environments.

CuroCell Cirrus® 2.0 combines air with a soft foam surface and provides good user comfort and pressure redistribution⁽²⁾. The added heel function offers pressure redistribution on the heels.

The silent and easy to use CuroCell Cirrus® pump with four different settings; alternating, constant low pressure, alternating + seat inflation and constant low pressure + seat inflation make this mattress functional. In addition the dual cover system reduces the shearing effect when adjusting the bed.

hygienic outer cover which conforms to relevant standards for hygiene, quality and safety⁽³⁾⁽⁵⁾. All covers are manufactured from modern, recyclable materials and are vapour permeable⁽⁴⁾, which contributes to reduce the risk of skin maceration.

REFERENCES

(1) European Pressure Ulcer Advisory Panel, National Pressure Injury Advisory Panel and Pan Pacific Pressure Injury Alliance. Prevention and Treatment of Pressure Ulcers/Injuries: Clinical Practice Guideline. The International Guideline. Emily Haesler (Ed). EPUAP/NPIAP/PPPIA: 2019.

(2) Examination of the shear force and the pressure relief capability, Berlin Cert 2009.

(3) ISO 16603 "resistance of penetration by blood and body fluids", ISO 16604 "resistance of penetration by blood-borne pathogens", EN 14126.

(4) SS-EN ISO 15496.

(5) SS 876 00 20:2017, Healthcare textiles - Mattresses - Specifications and requirements.



Always read the instruction for use prior to usage of the product.

Care of Sweden, CuroCell Cirrus are trademarks of Care of Sweden AB. © Care of Sweden AB, 2021. All rights reserved. Last update: 2021-06-07.

FEATURES & CHARACTERISTICS	FACTS
<ul style="list-style-type: none"> • Alternating pressure on soft foam surface • Cable holder • Digital and easy to use pump • Easy handling and cleaning • Heel slope, for pressure relief of heels • Hygienic cover with welded seams • Rapid deflation of the mattress • Silent and energy-saving pump • Soft, elastic, vapour permeable⁽⁴⁾ hygienic cover • Sound and light notification for low pressure • Sound and light notification for power failure • Caring mode (max firm) with automatic return • Transport function 	<ul style="list-style-type: none"> • Type: Dynamic mattress replacement system • Height: 21 cm • Function: an aid for prevention and treatment up to and including category 4⁽¹⁾ • Recommended user weight: up to 230 kg • Fire: EN 597-1, -2 • Power supply: 230 V/50 Hz • Noise level of pump: 27 dBA • Warranty: 2 years • Cleaning of cover: wipe with cleaning agent and/or disinfectants • Machine wash max 95 °C, tumble drying
SIZE	OPTIONAL
80/85/90/105/120x200 cm and 80/85/90x210 cm.	<ul style="list-style-type: none"> • Durable bag