

CuroCell® Nova CX13

Mattress overlay system with alternating pressure



CuroCell® Nova CX13 with cover Olivia

CuroCell® Nova CX13

The CuroCell® Nova CX13 is used as an aid in the prevention and treatment of pressure ulcers up to category 2⁽¹⁾. It is a mattress overlay system that is easy to handle and operate.

The CuroCell® Nova CX13 uses the alternating pressure method and is available in several different models. It has replaceable polyurethane cells, which are housed in a separate inner cover.

Also available is CuroCell® Nova CX19, a mattress replacement system equipped with an integrated foam core.

Care of Sweden's mattresses are supplied with a four way stretch hygienic outer cover which conforms to relevant standards for hygiene, quality and safety⁽²⁾⁽⁴⁾.



All covers are manufactured from modern, recyclable materials and are vapour permeable⁽³⁾, which contributes to reduce the risk of skin maceration.

REFERENSER

(1) European Pressure Ulcer Advisory Panel, National Pressure Injury Advisory Panel and Pan Pacific Pressure Injury Alliance. Prevention and Treatment of Pressure Ulcers/Injuries: Clinical Practice Guideline. The International Guideline. Emily Haesler (Ed). EPUAP/NPIAP/PPPIA: 2019

(2) CDC (2003), Guidelines for Environmental Infection Control in Health-Care facilities, updated 2019.

(3) SS-EN ISO 15496:2004, DIN 53122-1.

(4) SS 876 00 20:2017, Healthcare textiles - Mattresses - Specifications and requirements.

(5) EN ISO 11201 Acoustics - Noise emitted by machinery and equipment - Determination of emission sound pressure levels", SP 2018.



Always read the instructions for use prior to usage of the product.

Care of Sweden and CuroCell are trademarks of Care of Sweden AB.
© Care of Sweden AB, 2021. All rights reserved. Last update: 2021-06-03.

FEATURES & CHARACTERISTICS	FACTS
<ul style="list-style-type: none"> • CPR function – rapid deflation • Easy handling and cleaning • Replaceable cells in polyurethane • Silent, energy-saving pump • Soft, elastic, vapour permeable⁽³⁾ hygienic cover • Transport function 	<ul style="list-style-type: none"> • Type: alternating pressure mattress overlay • Height: 13 cm • Function: an aid for prevention and treatment up to and including category 2⁽¹⁾ • Recommended user weight: up to and including 150 kg • Input voltage: 230 V/50 Hz • Sound level, pump, max: 21 dBA⁽⁵⁾ • Warranty: 2 years • Machine wash max 95 °C, tumble drying • Cleaning of cover: wipe with cleaning agent and/or disinfectants
<p>OPTIONAL</p> <ul style="list-style-type: none"> • Durable bag • Heel function • Hygienic cover with welded seams 	
	<p>SIZES</p> <p>80/85/90/105/120x200 cm and 80/85/90x210 cm. If you need any other sizes please contact your distributor.</p>