

# CuroCell® Nova CX19

## Dynamic mattress replacement system



CuroCell® Nova CX19 with cover Olivia

FACTS	OPTIONAL
<ul style="list-style-type: none"> <li>• Type: Dynamic mattress replacement system</li> <li>• Height: 19 cm</li> <li>• Function: Prevention and treatment up to and including category 2<sup>(1)</sup></li> <li>• Recommended user weight: 0-150 kg</li> <li>• Fire: EN 597-1, -2</li> <li>• Input voltage: 230 V/50 Hz</li> <li>• Sound level, pump, max: 21 dBA</li> <li>• Warranty: 2 years</li> <li>• Machine wash max 95 °C, tumble drying</li> <li>• Cleaning of cover: Wipe with cleaning agent and/or disinfectants</li> </ul>	<ul style="list-style-type: none"> <li>• Durable bag</li> <li>• Sheet holder</li> <li>• Cable holder</li> <li>• Hygienic cover with welded seams</li> </ul> <p><b>SIZES</b></p> <p>80/85/90/105/120x200 cm and 80/85/90x210 cm. If you need any other sizes please contact your distributor.</p>

# CuroCell® Nova CX19

CuroCell® Nova CX19 is used as an aid in the prevention and treatment of pressure ulcers up to category 2<sup>(1)</sup>. It is used in all types of healthcare environments and is a mattress replacement system that is easy to handle and operate.


CuroCell® Nova CX19 works according to the alternating pressure principle and are available in several different versions. It has replaceable cells in polyurethane which are placed in a separate inner cover and is equipped with an integrated foam core.

The mattress is supplied with a four way stretch hygienic outer cover which conforms to relevant standards for hygiene, quality and safety<sup>(2)(4)</sup>. All covers are manufactured from modern, recyclable materials and are vapour permeable<sup>(3)</sup>, which contributes to reduce the risk of skin maceration.



## REFERENCES

- (1) European Pressure Ulcer Advisory Panel, National Pressure Injury Advisory Panel and Pan Pacific Pressure Injury Alliance. Prevention and Treatment of Pressure Ulcers/Injuries: Clinical Practice Guideline. The International Guideline. Emily Haesler (Ed). EPUAP/NPIAP/PPPIA: 2019.
- (2) CDC (2003), Guidelines for Environmental Infection Control in Health-Care facilities, updated 2019.
- (3) SS-EN ISO 15496:2004, DIN 53122-1.
- (4) SS 876 00 20:2017, Healthcare textiles - Mattresses - Specifications and requirements.

FEATURES CONTROL UNIT	FEATURES MATTRESS
<ul style="list-style-type: none"> <li>• Silent, energy-saving pump</li> <li>• Transport function</li> </ul>	<ul style="list-style-type: none"> <li>• CPR function – rapid deflation</li> <li>• Soft, elastic, vapour permeable<sup>(3)</sup> hygienic cover</li> <li>• Easy handling and cleaning</li> <li>• Heel function</li> <li>• Replaceable cells in polyurethane</li> </ul> <p> Always read the instruction for use prior to usage of the product.</p> <p>Care of Sweden and CuroCell are trademarks of Care of Sweden AB. © Care of Sweden AB, 2021. All rights reserved. Last update: 2021-06-07.</p>