

# Optimal Evacuation Sheet

For quick & simple evacuation

Optimal Evacuation Sheet makes it possible to quickly and simply move bedridden individuals away from danger. It can also be used for evacuation in stairs.

It is placed under the existing mattress and is fastened to the mattress with elastic straps. It is easy to embed and is available in several sizes to fit most mattresses.

Optimal Evacuation Sheet is ergonomically designed and produced from strong material that slides well on both wet and dry surfaces.

The Evacuation Sheet is X-ray transparent and equipped with handles for sideways movement.



Technical specification	
Material	Polyamide/polyester
Sizes (WxL)	80/90/105/120 cm x 200 cm
Fire	EN 597-1, -2
Standards	Fulfils requirements of the SS 8760039 (Swedish Standard) for rescue sheets.
Cleaning instructions	Wipe with cleaning agent and/or disinfectants. Machine wash max 95 °C, tumble drying.
CE-marking	Registered and marked in accordance with MDR (EU) 2017/745



Always read the instructions for use prior to use.

Care of Sweden and Curera are trademarks of Care of Sweden AB. © Care of Sweden AB, 2024. All rights reserved. Last update: 2024-01-22