

Optimal Hygiene Pillow

Comfortable & durable pillow

We offer a wide range of accessories, including pillows and pillow covers that are designed to meet healthcare needs and are also suitable for home environments.

Optimal Hygiene Pillow is a pillow that consists of ball fibers that can easily be fluffed up to maintain the original shape of the pillow.

Covers for comfort and hygiene

Optimal Hygiene Pillow comes with a liquidproof hygiene cover for easy cleaning.

The liquidproof hygiene cover is manufactured in a four-way stretch⁽¹⁾⁽²⁾ fabric to reduce the risk of shear forces, and is vapour-permeable⁽³⁾ to lower the risk of skin maceration. The hygiene cover also features a liquidproof zipper.



Technical specification	
Material	Optimal Hygiene Pillow Polyester/cotton Pillow Case Polyurethane/polyester
Sizes	35x45 cm 50x60 cm 50x70 cm
Cleaning instructions	Pillow Machine wash at max 60 °C Pillow case Wipe clean with detergent and/or disinfectant. Machine wash max. 95°C



References

- (1) CDC (2003), Guidelines for Environmental Infection Control in Health-Care facilities, updated 2019.
- (2) SS 876 00 20:2017, Healthcare textiles - Mattresses - Specifications and requirements.
- (3) SS-EN ISO 15496:2004, DIN 53122-1.



Always read the instructions for use prior to use.

Care of Sweden and Curera are trademarks of Care of Sweden AB. © Care of Sweden AB, 2026. All rights reserved. Last update: 2026-06-08.