

# PRIMA 6

A functional elevating pad



PRIMA 6 with cover Marin

# PRIMA 6

PRIMA 6 is an easy and functional basic mattress with good durability.

It is an elevating pad that is placed underneath a mattress overlay, for example CuroCell® Nova (air mattress) or Optimal 5zon® BM (preventative).

Care of Sweden's mattresses are supplied with a four way stretch hygienic outer cover which conforms to relevant standards for hygiene, quality and safety<sup>(1)(3)</sup>.

All covers are manufactured from modern, recyclable materials and are vapour permeable<sup>(2)</sup>, which reduce the risk of skin maceration.



## AVAILABLE COVERS

### Marin

With zipper on one side.

### Marin

With zipper on three sides.

## REFERENCES

(1) CDC (2003), Guidelines for Environmental Infection Control in Health-Care facilities, updated 2019

(2) SS-EN ISO 15496:2004, DIN 53122-1.

(3) SS 876 00 20:2017, Healthcare textiles - Mattresses - Specifications and requirements



Always read the instruction for use prior to usage of the product.

Care of Sweden, CuroCell and Optimal 5zon are trademarks of Care of Sweden AB. © Care of Sweden AB, 2023. All rights reserved. Last update: 2023-07-05.

FEATURES & CHARACTERISTICS	FACTS
<ul style="list-style-type: none"><li>• Easy to clean</li><li>• Hygienic cover</li><li>• Reversible</li></ul>	<ul style="list-style-type: none"><li>• Core: HR foam 35 kg/m<sup>3</sup></li><li>• Height: 6 cm</li><li>• Fire: EN 597-1, -2</li><li>• Warranty: 2 years</li><li>• Cleaning of cover: Wipe with cleaning agent and/or disinfectants</li><li>• Machine wash max 95 °C, tumble drying</li><li>• Prima 6 is designed to be used underneath our overlay mattresses. It can accommodate all patient weights applicable to our overlay mattresses.</li></ul>
SIZES	
80/85/90/105/120x200/210 cm. If you need any other sizes please contact your distributor.	