


# PRIMA BM

Mattress overlay that offers  
good comforty & durability



PRIMA BM with Cover Marin

# PRIMA BM

Mattress overlay PRIMA BM has a mattress core in one piece that provides the user good comfort and pressure relief. It is delivered with hygienic cover Marin.

Care of Sweden's mattresses are supplied with a four way stretch hygienic outer cover which conforms to relevant standards for hygiene, quality and safety<sup>(1)(3)</sup>.

All covers are manufactured from modern, recyclable materials and are vapour permeable<sup>(2)</sup>, which reduce the risk of skin maceration.



## AVAILABLE COVERS

### Marin

With zip on one side.

### Marin

With zip on three sides.

## REFERENCES

(1) CDC (2003), Guidelines for Environmental Infection Control in Health-Care facilities, updated 2019

(2) SS-EN ISO 15496:2004, DIN 53122-1.

(3) SS 876 00 20:2017, Healthcare textiles - Mattresses - Specifications and requirements



Always read the instruction for use prior to usage of the product.

Care of Sweden, CuroCell and Optimal 5zon are trademarks of Care of Sweden AB. © Care of Sweden AB, 2021. All rights reserved. Last update: 2021-05-28.

FEATURES & CHARACTERISTICS	FACTS
<ul style="list-style-type: none"><li>• Easy to clean</li><li>• Hygienic cover</li><li>• Reversible</li></ul>	<ul style="list-style-type: none"><li>• Core: HR foam 35 kg/m<sup>3</sup></li><li>• Height: 6 cm</li><li>• Recommended user weight: up to 130 kg*</li><li>• Fire: EN 597-1, -2</li><li>• Warranty: 2 years</li><li>• Cleaning of cover: wipe with cleaning agent and/or disinfectants.</li><li>• Machine wash max 95 °C, tumble drying.</li></ul> <p>*Please note that the underlying mattress affects the overall function including the recommended user weight.</p>
SIZES	
80/85/90/105/120x200 cm. If you need any other sizes please contact your distributor.	