

# Curera<sup>®</sup> Levelling plate

For levelling the seat surface in wheelchairs



# Curera<sup>®</sup> Levelling plate

If the seat surface in wheelchairs needs to be levelled, this is available as an accessory. This is a stable sheet which can be placed in the cover on the underside of the cushion. It is easy to clean and take care of.

Can be used with all Curera<sup>®</sup> seat cushions apart from Tremix which already has this function.



Always read the instructions for use prior to use.

Care of Sweden is a trademark of Care of Sweden AB. © Care of Sweden AB, 2023.  
All rights reserved. Last update: 2023-01-10.

AVAILABLE SIZES	FACTS
35x35 cm, 40x40 cm, 40x45 cm, 42x40 cm, 42x45 cm, 45x40 cm, 45x45 cm, 45x50 cm, 47x40 cm, 47x45 cm, 47x50 cm, 50x45 cm, 50x50 cm, 55x50 cm, 55x55 cm	<ul style="list-style-type: none"><li>• <b>Height: 2.5 cm (bevelled edges)</b></li><li>• <b>Warranty: 2 years</b></li><li>• <b>Materials: Cell Polyethylene</b></li><li>• <b>Classification: ISO-kod 181006-02</b></li><li>• <b>Fire: Complies with ISO 7176-16</b></li></ul>