

Optimal M3

Quality mattress replacement



Optimal M3 with Marin Cover.

FACTS	SIZES
<ul style="list-style-type: none">• Core: HR foam 50kg/m³• Height: 12 cm• Recommended user weight: Up to 130 kg• Fire: EN 597-1, -2• Warranty: 3 years• Cleaning of cover: Wipe with cleaning agent and/or disinfectants.• Machine wash max 95 °C, tumble drying.	<p>80/85/90/105/120x200 cm and 80/85/90x205/210 cm.</p> <p>If you need any other sizes please contact your distributor.</p>

Optimal M3

Optimal is a split core mattress with high density and is delivered with hygienic cover Marin and Marin/Evac.


The Marin cover is detachable with a zip fastener on one of its short sides. The Marin/Evac cover is equipped with a divisible, drip shield zip and feature handles for evacuation.

Care of Sweden's mattresses are supplied with a four way stretch hygienic outer cover which conforms to relevant standards for hygiene, quality and safety⁽¹⁾⁽²⁾. All covers are manufactured from modern, recyclable materials and are vapour permeable⁽³⁾, which contributes to reduce the risk of skin maceration.



REFERENCES

- (1) CDC (2003), Guidelines for Environmental Infection Control in Health-Care facilities, updated 2019
- (2) SS 876 00 20:2017, Healthcare textiles - Mattresses - Specifications and requirements
- (3) SS-EN ISO 15496:2004, DIN 53.122.

FEATURES & CHARACTERISTICS	AVAILABLE COVERS
<p>Pressure Elastic material with pressure equalisation that provides good comfort.</p> <p>Shear forces The mattress is specially designed to reduce the effects of shear forces.</p> <p>Moisture The hygienic cover is vapour permeable⁽³⁾, which contributes to reduce the risk of skin maceration.</p> <p>Hygiene The cover is easy to clean and has a good hygienic standard⁽¹⁾⁽²⁾.</p>	<p>Marin Drip shield zip on one side.</p> <p>Marin Drip shield zip on three sides.</p> <p>Marin Drip shield zip on three sides and handles on one side.</p> <p>Marin/Evac Divisible, drip shield zip and handles for evacuation on one side.</p> <p>In need of non hygienic cover, please contact your distributor.</p> <p> Always read the instruction for use prior to usage of the product.</p> <p>Care of Sweden and CuroCell are trademarks of Care of Sweden AB. © Care of Sweden AB, 2021. All rights reserved. Last update: 2021-06-09.</p>