

CuroCell® S.A.C

Seat cushion for support & comfort

CuroCell® S.A.C. is a static air seat cushion for comfort, prevention and as an aid in the treatment of pressure injuries up to category 2⁽¹⁾. It is pre-filled and easy to handle. The durable soft material for long term use provides a good comfort and support.



Always read the instructions for use prior to use.

REFERENCES

(1) National Pressure Injury Advisory Panel, European Pressure Ulcer Advisory Panel and Pan Pacific Pressure Injury Alliance. Prevention and Treatment of Pressure Ulcers/Injuries: Clinical Practice Guideline. The International Guideline: Fourth Edition. Emily Haesler (Ed.). 2025.

FEATURES AND CHARACTERISTICS

- Air venting holes, allows transportation of heat and moisture
- Durable soft material for long term use
- Light and easy to use
- Low profile design
- Patient ready
- Pre-inflated
- Reduces tissue interface pressure

FACTS

- Type: Static air seat cushion
- Function: Prevention and treatment up to and including category 2⁽¹⁾
- Recommended user weight: 0-160 kg
- Material: Film, color gray
- Warranty: 1 year
- Classification: ISO code 04 33 03 (pressure relieving)
- Sizes: 45x45x5 cm & 50x50x5 cm
- Cleaning: Wipe with cleaning agent and/or disinfectants

Care of Sweden is a trademark of Care of Sweden AB. © Care of Sweden AB, 2026.
All rights reserved. Last update: 2026-01-23.