

3 m-tubing kit

For easy access and increased independence

Enhancing independence and wellbeing

Supporting the independence of individuals with limited mobility can have a positive impact on both their physical and mental well-being. For those who can enter and exit the bed independently, having tools that simplify this process may reduce the need for assistance from healthcare staff. By making the control unit more accessible, it helps caregivers work more efficiently and allow users to manage their needs with greater ease, contributing to overall quality of care and independence.

When individuals can perform certain tasks independently, it frees up time for healthcare staff while giving the user more control and dignity in their daily life. For this reason, we have developed the 3 m-tubing kit, primarily designed to address the needs of the individual while also creating positive outcome for caregivers.

The 3 m-tubing kit is developed for individuals with limited mobility who can independently get in and out of bed. Compatible with the entire CX series, (including A4 and IQ control units). The tubing set offers greater flexibility, enabling the user to activate the MAX function for independent movement in and out of bed. Its extended length makes it easy to position the control unit on a bed rail or table, ensuring convenient access.

Secure and convenient with adjustable straps

3 m-tubing kit comes with two adjustable straps that make it easy to secure the tubing along the bed. By placing and fastening the straps at strategic points, the tubing is kept securely in place. This contributes to a safer and more functional environment for both the user and healthcare staff.

Technical specification	
Colour	Black
Size	3 m
Tubing protection material	100% Polyester Polyurethane
Compatible mattresses	CX10, CX15, CX16, CX20, CX20 Plus and Cirrus
Attachment solutions	2 adjustable straps for securing the tubing to the bed



Always read the instructions for use prior to use.

Care of Sweden and CuroCell are trademarks of Care of Sweden AB. © Care of Sweden AB, 2024. All rights reserved. Last update: 2024-11-21.