

Curera® Levelling plate

For levelling the seat surface in wheelchairs

If the seat surface in wheelchairs needs to be levelled, this is available as an accessory. This is a stable sheet which can be placed in the cover on the underside of the cushion. It is easy to clean and take care of.

Can be used with all Curera® seat cushions apart from Tremix which already has this function.



Always read the instructions for use prior to use.



Care of Sweden is a trademark of Care of Sweden AB. © Care of Sweden AB, 2025.
All rights reserved. Last update: 2025-05-08.

| AVAILABLE SIZES | FACTS |
|---|--|
| <p>35×35 cm, 40×40 cm, 40×45 cm, 42×40 cm, 42×45 cm, 45×40 cm, 45×45 cm, 45×50 cm, 47×40 cm, 47×45 cm, 47×50 cm, 50×45 cm, 50×50 cm, 55×50 cm, 55×55 cm</p> <p>More sizes can be ordered.</p> | <ul style="list-style-type: none">• Height: 2.5 cm (bevelled edges)• Warranty: 2 years• Materials: Cell Polyethylene• Classification: ISO-kod 122506-02• Fire: Complies with ISO 7176-16 |