

# Curera<sup>®</sup> Leg/Knee Separator

To separate your legs

The leg separator is used to keep your legs apart. It is made of sturdy foam with a degree of elasticity to provide support for the user.



Technical specification	
Sizes (WxH)	9/12/15 cm x 16 cm
Cleaning instructions	<p>Clean the leg separator with a mild detergent (e.g. washing-up liquid) and water or with a disinfectant/cleaner intended for this purpose such as DAX surface disinfectant or similar. Perform, Virkon or Desisoft can also be used.</p> <p>Let it dry in a warm, ventilated area (not in direct sunlight). It must be completely dry before it is used again.</p> <p>The cover is machine washable at 60°C and can be either tumble-dried or drip-dried.</p>
CE-marking	Registered and marked in accordance with MDR (EU) 2017/745



Always read the instructions for use prior to use.

Care of Sweden and Curera are trademarks of Care of Sweden AB. © Care of Sweden AB, 2026. All rights reserved. Last update: 2026-01-23.