

Optimal M3

Quality mattress replacement



Optimal M3 with Marin Cover.

FACTS	SIZES
• Core: HR foam 50kg/m³ • Height: 12 cm • Recommended user weight: Up to 180 kg • Fire: EN 597-1, -2 • Warranty: 3 years • Cleaning of cover: Wipe with cleaning agent and/or disinfectants. • Machine wash max 95 °C, tumble drying.	80/85/90/105/120x200 cm and 80/85/90x205/210 cm. If you need any other sizes please contact your distributor.

Optimal M3

Optimal is a split core mattress with high density and is delivered with hygienic cover Marin and Marin/Evac.

The Marin cover is detachable with a zip fastener on one of its short sides. The Marin/Evac cover is equipped with a divisible, drip shield zip and feature handles for evacuation.

Care of Sweden’s mattresses are supplied with a four way stretch hygienic outer cover which conforms to relevant standards for hygiene, quality and safety⁽¹⁾⁽²⁾. All covers are manufactured from modern, recyclable materials and are vapour permeable⁽³⁾, which contributes to reduce the risk of skin maceration.



REFERENCES

- (1) CDC (2003), Guidelines for Environmental Infection Control in Health-Care facilities, updated 2019
- (2) SS 876 00 20:2017, Healthcare textiles - Mattresses - Specifications and requirements
- (3) SS-EN ISO 15496:2004, DIN 53.122.

FEATURES & CHARACTERISTICS	AVAILABLE COVERS
Pressure Elastic material with pressure equalisation that provides good comfort. Shear forces The mattress is specially designed to reduce the effects of shear forces. Moisture The hygienic cover is vapour permeable⁽³⁾, which contributes to reduce the risk of skin maceration. Hygiene The cover is easy to clean and has a good hygienic standard⁽¹⁾⁽²⁾.	Marin Drip shield zip on one side. Marin Drip shield zip on three sides. Marin Drip shield zip on three sides and handles on one side. Marin/Evac Divisible, drip shield zip and handles for evacuation on one side. In need of non hygienic cover, please contact your distributor.  Always read the instruction for use prior to usage of the product. Care of Sweden and CuroCell are trademarks of Care of Sweden AB. © Care of Sweden AB, 2021. All rights reserved. Last update: 2026-10-01.